



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt Only for TDX® SP Power Wheelchair

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

Ordering a single power seating system on a multiple power option wheelchair base may require an affirmative determination through Advanced Determination of Medicare Coverage (ADMC). Please consult an attorney or other advisor to discuss this situation and any other specific billing situation in further detail. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program. This price list is for information purposes only and is not billing or legal advice.

POWERED SEATING BASE

TDXSP-CG TDXSP-CG for Formula CG with MK6i™ \$5,970.00 K0856

Group 3 - K0856 Basic Package	
TDXSP-CG TDXSP for Formula CG with MK6i	
300 lbs. Weight Capacity	Standard
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
4-Pole Motors	Standard
Flip Back Adjustable Height Arm	Standard
MK690ACC 90 Amp Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option CGT Tilt (E1002)	\$4,590.00

WHEEL OPTIONS

1430-3 14" x 3" Pneumatic Tires Standard
 1431-3 14" x 3" Tires with Foam Filled Inserts No Charge

FRAME FINISH

CONTEMPORARY

24P Wet Black No Charge
 60P Silver Metallic No Charge
 105P Midnight Blue No Charge
 115P Black Prism No Charge

TOO HOT

30P Sunny Yellow No Charge
 61P Electric Red No Charge
 119P Bubblegum Pink No Charge
 120P Grape Madness No Charge

WAY COOL

62P Electric Blue No Charge
 127P Grasshopper Green No Charge
 116P Sunkissed Orange No Charge
 117P Tangerine Red No Charge
 29P White Snow No Charge

OMIT SEAT FRAME

OSF Omit Seat Frame for Base Only Less \$400.00
 Base only to send for aftermarket seating.

TRANSPORT TIE DOWN

TRBKTS Wheelchair Transport Brackets \$250.00

TRANSPORT TIE-DOWN*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt Only for TDX® SP Power Wheelchair

SEAT

SYSTEM TYPE

- CGT CG Tilt Only System \$4,590.00 E1002
 CNVACT Conventional Actuator..... Standard

SEAT-TO-FLOOR HEIGHT

Nominal measurement at front edge of seat pan to floor on 18" deep seat.STF bracket allows for 1" of adjustability.

Due to variation in options the listed seat-to-floor heights have a tolerance of +/-0.5"

- LOW 17.25" to 18.25" at 0° No Charge
 MED 18.0" to 19.0" at 0° No Charge
 TALL 19.0" to 20.0" at 0° No Charge

SEAT SELECTION

- CTB Contoura™ Seating \$605.00 E2620 E2621
 Does not include seat cushion.
 BPO Back Posts Only/Round Canes Only No Charge
 CURV Curved No Charge

MATRX ELITE BACK 16"-20"w (3" contour)

- PBE14 Matrx Elite 14" high \$625.00 E2620
 PBE16 Matrx Elite 16" high \$625.00 E2620
 PBE 18 Matrx Elite 18" high \$625.00 E2620
 PBE 20 Matrx Elite 20" high \$625.00 E2620

SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range (Adjustable) No Charge
 2024W 20" to 22" Seat Width Range (Adjustable) No Charge

SEAT WIDTH SETTING

- W16 16" Seat WidthNo Charge
 W17 17" Seat Width (Not available with Contoura) No Charge
 W18 18" Seat WidthNo Charge
 W19 19" Seat Width (Not available with Contoura) No Charge
 W20 20" Seat WidthNo Charge
 W21 21" Seat Width (Not available with Contoura) No Charge
 W22 22" Seat Width No Charge

SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range (Adjustable) No Charge
 1922D 19" to 22" Seat Depth Range (Adjustable) No Charge

CONTOURA SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range (Adjustable) No Charge
 1922D 19" to 22" Seat Depth Range (Adjustable) No Charge

SEAT DEPTH SETTING

- D16 16" Seat DepthNo Charge
 D17 17" Seat DepthNo Charge
 D18 18" Seat DepthNo Charge
 D19 19" Seat DepthNo Charge
 D20 20" Seat DepthNo Charge
 D21 21" Seat Depth No Charge
 D22 22" Seat Depth No Charge

SEAT TILT

- 0SA 0°No Charge
 5SA 5°No Charge

ARM TYPE (PAIR)

- 19F Full Length Adj Hgt Conv. No Charge E0973
 29F Desk Length Adj Hgt Conv.No Charge E0973
 L8401 Omit ArmsNo Charge

FLIP UP MIXED ADJUSTABLE ARM TYPE (EACH)

All arms are detachable.

Right Arm

- 19FR Full Length Adj. Hgt. Right No Charge
 29FR Desk Length Adj. Hgt. Right No Charge

Left Arm

- 19FL Full Length Adj. Hgt. Left No Charge
 29FL Desk Length Adj. Hgt. Left No Charge

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS

- 7321BK Hook & Loop Chest Positioning Strap \$100.00 E0960

BACK HEIGHT

- TSBH20 Back Height 20" No Charge
 TSBH21 Back Height 21" No Charge
 TSBH22 Back Height 22" No Charge
 TSBH23 Back Height 23" No Charge
 TSBH24 Back Height 24" No Charge
 TSBH25 Back Height 25" No Charge
 TSBH26 Back Height 26" No Charge

BACK HEIGHT CONTOURA

- TSBH22 Back Height 22" No Charge
 TSBH25 Back Height 25" No Charge
 TSBH28 Back Height 28" No Charge

BACK ANGLE DEGREE

- 90BA 90° No Charge
 95BA 95° No Charge

Back Height Chart	
Back Type	Back Height
CTB - Contoura	22", 25" and 28"
BPO - Back Posts Only	16" - 24"
CURV - Curved	16" - 26"

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt Only for TDX® SP Power Wheelchair

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL)

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest \$210.00 E0955
- TAR127700HW Removable & Adjustable Hardware \$210.00 E1028

MATRX HEADREST PADS WITH STANDARD COVER

- MEHSP-6 Matrx Elan Standard Pad 6" \$235.00 E0955
- MEHSP-10 Matrx Elan Standard Pad 10" \$235.00 E0955
- MEHSP-14 Matrx Elan Standard Pad 14" \$235.00 E0955

MATRX HEADREST MOUNTING HARDWARE

- MEMAH Elan Multi-Axis Removable HDW \$220.00 E1028

LATERAL HIP SUPPORTS

HIP SUPPORT & PAD - LEFT

- LHSFL Hip Support with fixed HDW Left \$120.00 E0956
- LHS-XS Size XS (5.5"Wx4"H) No Charge
- LHS-S Size S (7"Wx6"H) No Charge
- LHS-M Size M (9"Wx4"H) No Charge
- LHS-L Size L (13"Wx4"H) No Charge

HIP SUPPORT & PAD - RIGHT

- LHSFR Hip Support with fixed HDW Right \$120.00 E0956
- RHS-XS Size XS (5.5"Wx4"H) No Charge
- RHS-S Size S (7"Wx6"H) No Charge
- RHS-M Size M (9"Wx4"H) No Charge
- RHS-L Size L (13"Wx4"H) No Charge

REMOVABLE HIP SUPPORT HARDWARE

Maxx Style Quick Release, Removable, Multi-Axis Mounting HDW for Hip Supports

- MQRHL MAXX Style Quick Release - Left \$215.00 E1028
- MQRHR MAXX Style Quick Release - Right \$215.00 E1028

ADDUCTORS

Hardware is packaged with equipment and cannot be ordered separately.

- ADD Adductors Small \$200.00 E0956
- ADDHW Hardware for ADD \$400.00 E1028

LATERALS (PAIR)

CONVENTIONAL BACK/BPO/CURVED BACK

Hardware is packaged with equipment and cannot be ordered separately.

- TCSLM Clamp-On Swingaway Laterals Medium \$210.00 E0956
 - TCSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TCSLL Clamp-On Swingaway Laterals Large \$210.00 E0956
 - TCSLLHW Removable and Adjustable Hardware \$210.00 E1028

LATERALS FOR CONTOURA SEATING

For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00
 - CSLMHW Removable and Adjustable Hardware \$210.00
- CSLL Swingaway Contoura Lat. Lg \$210.00
 - CSLLHW Removable and Adjustable Hardware \$210.00

MATRX FIXED LATERALS (FOR ELITE BACKS)

- PBLS-L Matrx Std Lateral - Fixed - Left (7.5"x5") \$100.00 E0956
- PBLS-R Matrx Std. Lateral - Fixed - Right (7.5"x5") .. \$100.00 E0956
- PBOL-L Matrx Offset Lateral Fixed Left \$110.00 E0956
- PBOL-R Matrx Offset Lateral Fixed Right \$110.00 E0956

MATRX SWING AWAY LATERAL PADS (FOR ELITE BACKS)

LEFT TRUNK SUPPORT PAD

- L-XS Size XS (3.25"Wx4.25"H) \$100.00 E0956
- L-S Size S (3.25"Wx5.25"H) \$100.00 E0956
- L-M Size M (3.75"Wx5.75"H) \$100.00 E0956
- L-L Size L (4.25"Wx6.5"H) \$100.00 E0956

RIGHT TRUNK SUPPORT PAD

- R-XS Size XS (3.25"Wx4.25"H) \$100.00 E0956
- R-S Size S (3.25"Wx5.25"H) \$100.00 E0956
- R-M Size M (3.75"Wx5.75"H) \$100.00 E0956
- R-L Size L (4.25"Wx6.5"H) \$100.00 E0956

SWING AWAY HARDWARE (FOR ELITE BACKS)

- SAEHBL Matrx Elite Swing Away HDW Left \$195.00 E1028
- SAHEBR Matrx Elite Swing Away HDW Right \$195.00 E1028
- OSAHL Matrx Offset Elite Swing Away HDW Left ... \$195.00 E1028
- OSAHR Matrx Offset Elite Swing Away HDW Right . \$195.00 E1028

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

ELECTRONICS

MK6i™ CONTROLLER

- MK690ACC Controller with Actuator Control Standard
Select this controller for all single or 2 actuator systems except for systems with recline and power center mount.
- MK690 90 AMP Controller No Charge
Select this controller for all systems with 3 or more actuators and with recline and power center mount systems.
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides Improved Tracking Capabilities and Enhanced Traction. Does not require Special Motors.

SINGLE ACTUATOR SYSTEMS - NON EXPANDABLE

- SPJAP SPJ™ + for Power Seating Standard
Includes Egg switch with 611 mount to operate single actuator.

ACTUATOR OPERATION THROUGH SPJ+

- SPJAPM SPJ+ with Joystick Actuator Control \$1,300.00 E2310
Allows actuator operation through the Driver Control. Includes port for the Attendant Switch.

EXPANDABLE CONTROLLER SYSTEM

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+ Multiple Drive Joystick with Color Display .. \$750.00 E2377
- ASLRDYM6 MK6i Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.

FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325

SIP-N-PUFF KIT

- PKG32666 Therafin Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- PACM6 Prop Attendant Control with quick release mount \$850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 E1028
QLAM6 includes ARM250.
- QLA Quad Link Retract. Joystick Mnt SPJ+ \$279.00 E1028
QLAM6 includes ARM250.
- ARM250 Stealth Height Adjustable Joystick Mount \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USB Ready) \$90.00
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 and 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 and 4 \$795.00
- MPY Monoport "Y" Cable for CMPJ+ or Display \$15.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt Only for TDX® SP Power Wheelchair

POWERED SEATING ACCESSORIES

- SANODE One Act. Oper. through Driver Ctrl \$1,300.00 E2310
Required for Tilt Operation through Driver Control - Single Actuator Systems only. Not available with SPJ+ Joystick. Does not come with Egg switch.
- 4WSB Four Way Switch Box No Charge
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge

Systems

For systems with 2 actuators NOT operating through driver control, choose 4WSB with either 4-way toggle or 4-push button.

For systems with 2 actuators OPERATING THROUGH DRIVER CONTROL, choose S4WSB.
Uses 4-way toggle or 4-push button as attendant control.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black \$156.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS

BATTERY TRAY TYPE

22TRY Base Width 24": 24TRY Base Width 25 1/2"

- 22TRY 22NF Style Tray Standard
- 24TRY 24 Group Style Tray No Charge

BATTERIES

- 22NBATTERY 22NF Battery (On Chair) \$300.00 E2361
- 24BATTERY 24Gel Battery (On Chair) \$410.00 E2363

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger No Charge
- L8900 Omit Battery Charger Less \$80.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00

FRONT RIGGINGS

HANGER BRACKET TYPE

All front riggings include a standard footplate.

- ASB Center Pivot Style Bracket No Charge

CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests 13.5" - 17.5" Standard
- 70STAPER 70° Tapered Footrests 13" - 17" No Charge
Standard with impact guards.
- AT5044 Manual Elevating Legrests \$255.00 E0990

CENTER MOUNT

- CTMT-F Center Mount Adj. Knee Angle (12.25" - 15") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.
- SML Small Footplates (5.25" x 7.5") No Charge
- MDM Medium Footplates (5.5" x 9.5") No Charge
- LRG Large Footplates (5.90" x 11.5") No Charge
- CALF Calf Pads For Center Mount \$200.00
- CTMT-FJ Jr. Center Mount Adj. Knee Angle (9.25" - 12") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.
- SML Small Footplates (5.25" x 7.5") No Charge
- MDM Medium Footplates (5.5" x 9.5") No Charge
- LRG Large Footplates (5.90" x 11.5") No Charge
- CALF Calf Pads For Center Mount \$200.00

POWER FOOT POSITIONING OPTIONS

LNX and LNX-L come Standard with curved individual calf pads that are independently height and angle adjustable. KtH range: LNX=13" to 15.5", LNX-L=13" to 19"

- LNX Power Center Mount Foot Platform \$2,995.00 E1012
Comes standard with a footboard. Not available with SPJ+ Joystick. Clear
- SML Small Footplates (5.25" x 7.5") \$275.00
SML comes with 2" of width adjustment ranging from 12-14"
- MDM Medium Footplates (5.5" x 9.5") \$275.00
MDM comes with 2" of width adjustment ranging from 13-15"
- LRG Large Footplates (5.9" x 11.5") \$275.00
LRG comes with 2" of width adjustment ranging from 14-16"
- LNX-L Power Center Mount Foot Platform- Long \$2,995.00 E1012
If 0SA is selected, LNX-L is only available with Tall STF height. Comes standard with a footboard. Not available with SPJ+ Joystick.
- SML Small Footplates (5.25" x 7.5") \$275.00
SML comes with 2" of width adjustment ranging from 12-14"
- MDM Medium Footplates (5.5" x 9.5") \$275.00
MDM comes with 2" of width adjustment ranging from 13-15"
- LRG Large Footplates (5.9" x 11.5") \$275.00
LRG comes with 2" of width adjustment ranging from 14-16"

FOOTPLATE OPTIONS (PAIR)

- 1651 Flip - Up Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$51.00
- HLOOP Heel Loops (pair) \$45.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula™ CG Tilt Only for TDX® SP Power Wheelchair

Tech Tips

Operating Powered Seating through the Driver's Control.
All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) - E2310.

Single Actuator Systems with CMPJ+, MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310.

Multiple Actuator Systems:

Use of the M6Ki enables actuator operation through the driver control.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Batteries (not included)

22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) MUST be used. GP24 batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations MUST be followed otherwise injury and damage may occur.

Contoura Seat Sizes

	D16	D17	D18	D19	D20	D21	D22
W16	*	*	*	*	*	*	*
W18	*	*	*	*	*	*	*
W20	*	*	*	*	*	*	*
W22	*	*	*	*	*	*	*

*Not available with Contoura Seat Cushion

Back Heights 22", 25" or 28"

Curved or BPO Seat Sizes

Seat Widths	16" - 22"
Seat Depths	16" - 22"
Back Heights	16" - 26"

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Formula CG Tilt Only Features and Specifications

- Weight capacity of 300 lbs.
- Adjustable seat frame

Comments

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way

Elyria, OH USA

44035-4190

800-333-6900

©2018 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.



Yes, you can.®

Revised July 2, 2018