



# Freedom<sup>®</sup> NXT<sup>™</sup>

Generation NEXT **Mini** Tilt-in-Space Mobility

Freedom: for the littlest Kids  
(and their best friends too)



NEW height adjustable  
stroller push handles

Super  
lightweight  
yet strong  
& dependable

Easy to adjust  
for growing kids

Frame widths as narrow as 10" and frame depths as short as 8"

Good things do come in small packages!

**NXT<sup>™</sup>Mini**  
Generation NEXT



# What makes the Mini so mighty...

## For little kids?

The NXT Mini is designed to fit the "little" kids more precisely. The seat depth adjusts 4" in ½" increments from only 8" to 12." With Freedom Seating designed for "Toddlers", the Mini can accommodate as little as a 5" seat depth. For accurate foot positioning, footrest heights can start at 3" from the top of the seat frame and adjust in ½" increments to 10" with a 4" extension. The new flip back armrest is only 8.5" long to support little arms and we now have smaller foot plates for supporting little feet.



The Mini is designed with little people's safety in mind

## For Moms and Dads?

The NXT Mini is packed with features to make using a wheelchair safe and simple. Freedom Designs has the only attendant hand controlled combination tilt and braking system. This system controls the speed on un-level surfaces, helps turn the wheelchair, position the seat frame in and out of tilt and locks the wheels, all from the attendant's location. Height adjustable stroller push handles fit the tallest moms and the shortest dads and when the frame is folded for transport, they fold down flat so the Mini fits into tiny little spaces. The anti-tippers are spring activated so that when pushing the wheelchair, they can be rotated forward under the frame and out of the way. For families on the go, the Mini folds into a super lightweight compact package that is easy to fold, lift and stow for travel.

## NXT Mini Specifications

### Seat Frame Width:

Adjusts in width without parts from 10" to 12" in 1" increments

### Seat Frame Depth:

Adjusts in depth without parts from 8" to 12" in ½" increments

### Seat Frame Heights:

15" to 19" with various combinations of casters and wheels

### Foot Plate Heights:

3" to 10" in ½" increments with various mounts or extensions

### Height Adjustable Stroller /Push Handle:

Standard is 18" to 24" and tall is 20" to 26"

### Tilt Range:

0° to 45° with ability to set max tilt angle at anything in between

### Tilt Operations:

Dual hand lever, foot activated or combination hand activated tilt, brake and wheel lock system

### Brake Operations:

Standard manual, foot activated, hand activated or combination hand activated tilt, brake and wheel lock system

### Rear Wheels:

12", 16" & 18" Mag and 16" & 18" Spoke

### Front Casters:

3", 4", 5", 6" & 7" in various configurations

### Rider Weight Capacity:

75 lbs.

### Frame Transport Weight:

19 lbs. in its folded transport configuration

### Frame Colors:

25 fun colors with seating upholstery to match



Our kids are our passion.™

**Freedom Designs Incorporated**  
2241 N. Madera Road, Simi Valley, CA 93065

phone: 800-331-8551

fax: 888-582-1509

e: CustomerService@freedomdesigns.com

www.freedomdesigns.com

Freedom Designs Incorporated (FDI) maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without notification. Check with FDI for the latest information.

Copyright © 2014 Freedom Designs, Inc.  
All rights reserved.  
Printed in the USA.  
WBNXM 112014



Manufactured in the U.S.A.