



Contour Static Seating on ROVI X3 Base
US PRICE LIST AND ORDER FORM
 Price Effective February 1, 2019

Reset Form




REQUIRED INFORMATION

Date of Order: _____ Dealer Account #: _____
 Dealer Name: _____
 PO #: _____ Tag _____
 After Market Items shipping to Motion Concepts: Yes No
 List Items: _____
 PO #: _____
 Name of the Manufacturer: _____
 Purchasing Contact: _____
 Phone: _____ Fax: _____
 E-mail: _____
 ATP/Therapist: _____
 Ship to Address: _____
 City: _____ State: _____
 Zip: _____ Client Gender: M F
 Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	



Weight Capacity is up to 300lbs unless specified otherwise



HCPCS CODES	K0848/, K0849
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH	23.25"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"
CASTER WHEEL SIZE (front/rear)	6"
TURNING RADIUS	19.5"
DRIVETRAIN	4-pole EAD Motors

BATTERIES	2 x M34
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs (M34)
RANGE ⁽²⁾	17.7 miles
WEIGHT CAPACITY	300 lbs
BASE WEIGHT (w/o batteries)	157 lbs
GROUND CLEARANCE	3"
MAXIMUM SPEED	6.2 mph

BASE VITALS

⁽¹⁾ All X3 wheels come with solid inserts

⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SOLID SEAT PAN - K0848



Item # ROVI X3-SS or Item # 888WNLEX3-SP-SS

CAPTAIN SEAT - K0849



Item # ROVI X3-CS or Item # 888WNLEX3-SP-CS

Contoured Static Seating on ROVI X₃ Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

BASE COLOR

Select Color: Standard With Lights⁽¹⁾ **\$1,195**

- Blizzard White Part # 800914-00200-03 Part # 800914-00201-03
- Matte Black Part # 800914-00200-50 Part # 800914-00201-50
- Solar Yellow Part # 800914-00200-69 Part # 800914-00201-69
- Motion Green Part # 800914-00200-36 Part # 800914-00201-36
- Original Orange Part # 800914-00200-23 Part # 800914-00201-23
- Bright Blue Part # 800914-00200-40 Part # 800914-00201-40
- Red Delicious Part # 800914-00200-21 Part # 800914-00201-21
- Party Pink Part # 800914-00200-2C Part # 800914-00201-2C
- Purple Reign Part # 800914-00200-80 Part # 800914-00201-80

Matching Battery Doors (select fender color first) . . . **\$495**

NOTES:

1. Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. Must select Joystick for Lights feature or Omni for operation.



Please see Custom Design Options on page 5/6

TIRE OPTIONS

- Standard Solid Tires Part # STIRES **STD**
- Pneumatic Tires Part # PTIRES **\$199**

MILITARY DECAL

3" Diameter Decal (Uninstalled) **\$50**



- Army Navy
- Air Force Marines

! Please note that the HCPCS codes indicated on this order form are our interpretation of the new codes and how we believe our products should be coded. The individual DMERCs will determine if our interpretation is appropriate and whether or not items coded will be accepted or denied. Motion Concepts cannot guarantee that use of the designated codes will ensure coverage or payment. Contact the insurer to confirm compliance with coding and coverage policy

FULL CAPTAIN-STYLE SEAT - K0849

All bases below come standard with Full Captain-Style Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

ROVI X₃⁽¹⁾ HCPCS Code: K0849 Part # 888WNLEX3-SP-CS. **\$6,695**
Select One:

- Standard 5.0 mph High Torque Motors** **N/C**
Part # HTM5
- Optional 5.8 mph Quick Speed Motors⁽²⁾** **\$200**
Part # QSM58
- Optional 6.2 mph High Speed Motors⁽²⁾** **\$1,550**
with OnTraxx Enhanced Driving Module Part # HSM62

- Seat Size: Grey** 18"x18" Part # SV18171/222RG
Black 18"x18" Part # SV18171/222RB

STF for 18 x 18⁽²⁾: 19" 20" 21" 22"

- Seat Size: Grey** 20"x20" Part # SV2121211/2RG
Black 20"x20" Part # SV2121211/2RB

STF for 20 x 20⁽²⁾: 20" 21" 22" 23"

NOTES:

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt, Weight capacity is 300 lbs.
2. Seat to floor height is measured at the seam of cushion on 20" widths and at the top of the leg trough on 18" widths. Seat to floor heights may vary +/- .25" dependent on base conditions.

SOLID SEAT PAN & CAPTAIN-STYLE SEAT - K0848

All bases above come standard with a Solid Seat Pan style Captain Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

ROVI X₃⁽¹⁾ HCPCS Code: K0848 Part # 888WNLEX3-SP-SS. **\$6,550**
Select One:

- Standard 5.0 mph High Torque Motors** **N/C**
Part # HTM5
- Optional 5.8 mph Quick Speed Motors⁽²⁾** **\$200**
Part # QSM58
- Optional 6.2 mph High Speed Motors⁽²⁾** **\$1,550**
with OnTraxx Enhanced Driving Module Part # HSM62

- Select:** **Grey** Seat Size: 18"x18" Part # SP181821RG
 20"x20" Part # SP202023RG
 Black Seat Size: 18"x18" Part # SP181821RB
 20"x20" Part # SP202023RB

Select STF⁽²⁾: 17" 18" 19" 20"


NOTES:

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt, Weight capacity is 300 lbs.
2. Seat to floor height is measured at the front of the seat pan based on a standard seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18".

Contoured Static Seating on ROVI X3 Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

SEAT CUSHION REFERENCE GUIDE



Cushion	Meshtex	Startex	Infection	Spacetex	Fusion	Onyx
Libra	N/A	STD	\$50	N/C	N/A	N/A
PS	N/A	STD - Reversible	STD - Reversible	N/A	N/A	N/A
PSP	N/A	STD - Reversible	STD - Reversible	N/A	N/A	N/A
Vi	STD	\$50	\$50	N/A	\$150	\$100

SEAT CUSHION OPTIONS - FOR SOLID SEAT

- Matrx Libra Seat Cushion 18" w-20" w** HCPCS Code E2624 . **\$485**
- Matrx PS Seat Cushion 18" w - 20" w** HCPCS Code E2605 . **\$345**
- Matrx PSP Seat Cushion 18" w - 20" w** HCPCS Code E2607 . **\$395**
- Matrx Vi Seat Cushion 18" w - 20" w** HCPCS Code E2607 . **\$425**

FABRIC OPTIONS FOR CUSHIONS

- Spacetex Fabric Option for Libra Cushion Part # SFLF N/C
- Infection Control Fabric Option for Libra Cushion Part # 10FL \$50
- Fusion Fabric Upgrade for Matrx Vi Cushion Part # FFV . . . \$150
- Onyx Fabric Upgrade for Matrx Vi Cushion Part # 0CUMA . . . \$100
- Startex Fabric Upgrade for Matrx Vi Cushion
- Smooth Side up Part # SFVS \$50
- Fabric Side up Part # SFVF \$50

SINGLE FUNCTION CONTROL OPTIONS

- VR2 Non-Expandable Joystick** Part # D50680 **STD**
- R-Net LED Remote Joystick⁽¹⁾⁽²⁾** Part # D51122SC **\$295**
- R-Net LED Remote Joystick for Lights⁽¹⁾** Part # D51108 **\$500**
- R-Net Color Remote Joystick with Monojacks⁽¹⁾⁽²⁾** Part # D51165 **\$600**
- R-Net Color Remote Joystick with Monojacks for Lights⁽¹⁾** Part # D50982 **\$650**
- R-Net Advanced Joystick⁽¹⁾** Part # D51711 **\$1295**
3.5" high-resolution VGA color display. IR Control as standard Includes Integrated Bluetooth module.
- R-Net Expandable Controller⁽³⁾** **\$615**
HCPCS Code: E2377 Part # D51109
- R-Net Harness for Expandable Controller⁽³⁾** **\$415**
HCPCS Code: E2313 Part # REH

- NOTES:**
1. Requires an R-Net Controller and Harness
 2. Light package not available with the joystick.
 3. Required for R-Net Remote Joysticks.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

JOYSTICK MOUNT

Fixed for Solid Pan/Captian Seat Part # 1343-00707-00 . . . **N/C**
 Select Mounting Position: Left Right
 Part # ML Part # MR

Swing Away Mount Bracket **\$295**
 Select Mounting Position: Left Part # 1343-00709-00-A3
 Right Part # SAQR

JOYSTICK HANDLES

- "T" Handle for Joystick Knob** **\$115**
HCPCS Code: E2323 Part # SA80941
- Ball Handle for Joystick Knob** **\$115**
HCPCS Code: E2323 Part # SA72441
- BodyPoint Joystick Handle** HCPCS Code: E2323 **\$120**
 - 3" U-Shaped Part # PC101B-20
 - 4" U-Shaped Part # PC102B-20
 - Domed Rubber Part # PC107B-20

ELECTRONIC OPTIONS

- R-Net Output Module(EMC)⁽¹⁾⁽²⁾** **\$995**
Part # D50883
- R-Net Bluetooth Module⁽¹⁾⁽²⁾** **\$895**
Select OS: iOS Apple Part # D51580 Android/PC Part # D51111
- Smart Phone Holder** Part # SPH **N/C**

- NOTES:**
1. Requires an R-Net Controller and Harness order above
 2. Must order R-Net Color Remote joystick with these options

ARMRESTS (MUST BE SELECTED)

- Adjustable Height Flip Back Arm Left** Part # 104300-00112. . **\$125**
- Adjustable Height Flip Back Arm Right** Part # 104400-00112 **\$125**

OTHER OPTIONS

- Wheelchair Transport Brackets⁽¹⁾** Part # XBKTS **\$250**
- NOTES:**
1. For unoccupied use only

BATTERIES

- 2 x MK GROUP 34 Batteries** Part # XBKTS M34 SLD G(Two) **\$650**
HCPCS Code: E2359

Total Retail Value \$ _____

Email Form

Contoured Static Seating on ROVI X3 Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CUSTOM FENDER DESIGN BY BOB RICE

Custom Fender Design by Bob Rice \$500

Standard Fender Configuration Part # Below N/C

Fenders w/Lights Configuration Part # Below + "L" . . . \$1195



Flames
Part # CFF



Personalized
Part # CFP



Camouflage
Part # CFC



Silver Concrete
Part # CFSC



U.S. Navy
Part # CFUSN



U.S. Coast Guard
Part # CFUSCG



U.S. Air Force
Part # CFUSAF



U.S. Marines
Part # CFUSAM



U.S. Army
Part # CFUSAA

DESIGN DESCRIPTION

(PLEASE USE IF YOU SELECTED PERSONALIZED OPTION)

Large empty box for design description.

DEALER SIGN OFF

This concept is agreed to and accepted by:

Dealer Name: _____

PO #: _____

Approved with no changes

Approved with changes

Review again

Date: _____

Signature: _____

Contoured Static Seating on ROVI X3 Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CUSTOM DIPPED DESIGN OPTIONS

Please select Fender option first, if battery doors are to match the fenders. If not make the design selection below.

Dipped Fender Designs

- Standard Fender Configuration**
Part # CDF + Code..... **\$395**
- Fenders w/Lights Configuration**
Part # CDFL + Code..... **\$1,590**
- Front Battery Door** Part # FBD + Code..... **\$395**
- Rear Battery Door** Part # RBD + Code..... **\$395**
- Front/Rear Battery Doors**
Combo Credit Part # FRBDC..... **(\$40)**
- Fenders & Front/Rear Battery Doors**
Combo Credit Part # FFRBDC..... **(\$85)**



CUSTOM DIPPED DESIGN OPTIONS



- U.S. Flag**
Code # 1
- Camouflage**
Code # 2
- Electric Circuit**
Code # 3
- Flames**
Code # 4
- 80's Green Rock**
Code # 5
- 80's Black Rock**
Code # 6
- Pink Camo**
Code # 7
- Flower Power**
Code # 8
- Graffiti**
Code # 9
- Animal Paw Print**
Code # 10