



**Invacare® ProSPIN™ X4™ Ultralightweight  
Active Performance Folder  
AADL Customers Only  
Effective October 1st, 2019**

NAME \_\_\_\_\_ PHONE NO \_\_\_\_\_  
 ACCOUNT # \_\_\_\_\_ P.O. \_\_\_\_\_ DATE \_\_\_\_\_  
 ADDRESS \_\_\_\_\_

**FRAME**

**FRAME TYPE**

- PROX4S X4 Folding Frame w/Swingaway Front End \$2,569
- PROX4F70 70 Degree Fixed Front Frame \$3,009  
*Not available with Heavy Duty, TRRO or TRBKTS.*
- PROX4F80 80 Degree Fixed Front Frame \$3,009  
*Not available with Heavy Duty, TRRO or TRBKTS.*

**FRAME SELECTIONS**

**FRAME WIDTH**

- W14 14" Wide See Matrix  
*Not available with HD - Heavy Duty Option*
- W15 15" Wide See Matrix  
*Not available with HD - Heavy Duty Option*
- W16 16" Wide See Matrix  
*Not available with HD - Heavy Duty Option*
- W17 17" Wide See Matrix  
*Not available with HD - Heavy Duty Option*
- W18 18" Wide See Matrix
- W19 19" Wide See Matrix
- W20 20" Wide See Matrix
- W22 22" Wide (22"D only) \$139  
*Must order HD - Heavy Duty Option.*

**SEAT DEPTH**

- D14 14" Deep See Matrix  
*Not available with HD - Heavy Duty Option*
- D15 15" Deep See Matrix  
*Not available with HD - Heavy Duty Option*
- D16 16" Deep See Matrix  
*Not available with HD - Heavy Duty Option*
- D17 17" Deep See Matrix  
*Not available with HD - Heavy Duty Option*
- D18 18" Deep See Matrix
- D19 19" Deep See Matrix
- D20 20" Deep (18-20"W only) \$139
- D21 21" Deep (20"W only) \$139  
*Must order HD - Heavy Duty Option.*
- D22 22" Deep (20" & 22"W only) see Matrix for other widths \$139  
*Must order HD - Heavy Duty Option.*

**HEAVY DUTY OPTION**

- HD Heavy Duty Option \$600  
*Only available on the PROX4S, X4 Folding Frame w/Swingaway Front End. Not Available with TRBKTS or TRRO option.*

**FRONT SEAT TO FLOOR**

- FSTF155 Seat-to-floor 15.5"
- FSTF16 Seat-to-floor 16"
- FSTF165 Seat-to-floor 16.5"
- FSTF17 Seat-to-floor 17"
- FSTF175 Seat-to-floor 17.5"
- FSTF18 Seat-to-floor 18"
- FSTF185 Seat-to-floor 18.5"
- FSTF19 Seat-to-floor 19"
- FSTF195 Seat-to-floor 19.5"
- FSTF20 Seat-to-floor 20"
- FSTF205 Seat-to-floor 20.5"
- FSTF21 Seat-to-floor 21"
- FSTF215 Seat-to-floor 21.5"

RSTFH/FSTFH COMPATIBILITY																
RSTFH	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
FSTFH																
15.5"	Y	Y	Y	Y	Y	N	N	N	N	N	N	N	N	N	N	N
16"	Y	Y	Y	Y	Y	Y	N	N	N	N	N	N	N	N	N	N
16.5"	N	Y	Y	Y	Y	Y	N	N	N	N	N	N	N	N	N	N
17"	N	N	Y	Y	Y	Y	Y	N	N	N	N	N	N	N	N	N
17.5"	N	N	N	Y	Y	Y	Y	Y	Y	N	N	N	N	N	N	N
18"	N	N	N	N	Y	Y	Y	Y	Y	N	N	N	N	N	N	N
18.5"	N	N	N	N	N	Y	Y	Y	Y	Y	N	N	N	N	N	N
19"	N	N	N	N	N	N	Y	Y	Y	Y	Y	Y	N	N	N	N
19.5"	N	N	N	N	N	N	N	Y	Y	Y	Y	Y	Y	N	N	N
20"	N	N	N	N	N	N	N	N	Y	Y	Y	Y	Y	Y	Y	N
20.5"	N	N	N	N	N	N	N	N	N	Y	Y	Y	Y	Y	Y	Y
21"	N	N	N	N	N	N	N	N	N	N	Y	Y	Y	Y	Y	Y
21.5"	N	N	N	N	N	N	N	N	N	N	N	Y	Y	Y	Y	Y

**REAR SEAT TO FLOOR**

- RSTF135 Seat-to-floor 13.5" No Charge
- RSTF14 Seat-to-floor 14" No Charge
- RSTF145 Seat-to-floor 14.5" No Charge
- RSTF15 Seat-to-floor 15" No Charge
- RSTF155 Seat-to-floor 15.5" No Charge
- RSTF16 Seat-to-floor 16" No Charge
- RSTF165 Seat-to-floor 16.5" No Charge
- RSTF17 Seat-to-floor 17" No Charge

**\* TRANSPORT TIE-DOWN**

TRRO includes four factory-installed transport brackets and a wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI/RESNA WC Vol 1. Section 19 Front Impact Test requirements for Wheelchairs with a 168 lb. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs.

TRBKTS includes four factory-installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

**\* Transport Tie-Down option is not retrofittable to existing models and is not field serviceable.**

Width & Depth Matrix									
Depth	14"	15"	16"	17"	18"	19"	20"	21"	22"
Width									
14"	X	X	X	X	X				
15"	X	X	X	X	X				
16"	X	X	X	X	X				
17"	X	X	X	X	X				
18"	X	X	X	X	X*	X*	X*		
19"			X	X	X*	X*	X*		
20"			X	X	X*	X*	X*	*	*
22"									*

**X - Available Sizes**  
**X\* - Heavy Duty Optional**  
**\* - Heavy Duty Only**

**REAR SEAT TO FLOOR - CONTINUED**

- RSTF175 Seat-to-floor 17.5" No Charge
- RSTF18 Seat-to-floor 18" No Charge
- RSTF185 Seat-to-floor 18.5" No Charge
- RSTF19 Seat-to-floor 19" No Charge
- FSTF195 Seat-to-floor 19.5" No Charge
- RSTF20 Seat-to-floor 20" No Charge
- FSTF205 Seat-to-floor 20.5" No Charge

**SUSPENSION SELECTIONS**

**TRANSPORT TIE DOWN**

- TRRO Transport Ready Option \$416
- TRBKTS Wheelchair Transport Brackets No Charge

**FRONT SUSPENSION (SOFR add 1" to selected FSTF & RSTF)**

- SOFR-F Frog Legs® Front Suspension Forks \$428

**FRONT SUSPENSION SETTING (MUST select SOFR-F above to pick your setting)**

- SET160 Soft Setting (Max. User Weight 160 lbs) No Charge
- SET300 Firm Setting (Max. User Weight 300 lbs) **Standard**

**ARM SELECTIONS**

**ARM TYPE**

- AT99 Swingaway Padded Arm \$261  
*Not available with Heavy Duty, Angle Adj. (AAB126 & AAB1720), TRRO, TRBKTS or ADJTENS.*
- CAD Conventional Adjustable Ht Desk Lgth Arms (pair) \$372
- CAF Conventional Adjustable Ht Full Lgth Arms (pair) \$372
- TAF T-Arms, Full, Adjustable Height \$372  
*Rigid sideguards included with the T-Arms. Bolt-on attachment*
- TAD T-Arms, Desk, Adjustable Height \$372  
*Rigid sideguards included with the T-Arms. Bolt-on attachment*
- LESSARMS Less Arms No Charge

**T-ARM HEIGHT**

- L6-10 Low Range 6" - 10" (adjusts in 1/2" inc) No Charge
- H8-14 High Range 8" - 14" (adjusts in 1" inc) No Charge

**BACK SELECTIONS**

**BACK TYPE**

- BACK1014 Adj. Ht. Back w/o Push Handles (10" - 14") No Charge
- BACK1216 Adj. Ht. Back w/o Push Handles (12" - 16") No Charge
- BACK1418 Adj. Ht. Back w/o Push Handles (14" - 18") No Charge
- BACK1620 Adj. Ht. Back w/o Push Handles (16" - 20") No Charge
- PH1216 Adj. Ht. Back w/ Push Handles (12" - 16") \$61
- PH1418 Adj. Ht. Back w/ Push Handles (14" - 18") \$61
- PH1620 Adj. Ht. Back w/ Push Handles (16" - 20") \$61
- A1418 10° Bend Back w/ Push Handles (14"-18") No Charge
- A1620 10° Bend Back w/ Push Handles (16"-20") No Charge
- AAB1216 Adjustable Angle Back (12"-16") \$172
- AAB1720 Adjustable Angle Back (17"-20") \$172

**BACK HEIGHT**

- ABH10 10" High No Charge
- ABH11 11" High No Charge
- ABH12 12" High No Charge
- ABH13 13" High No Charge
- ABH14 14" High No Charge
- ABH15 15" High No Charge
- ABH16 16" High No Charge
- ABH17 17" High No Charge
- ABH18 18" High No Charge
- ABH19 19" High No Charge
- ABH20 20" High No Charge

**BACK UPHOLSTERY TYPE**

- SLIPBACK Slip On Back Upholstery (Black) **Standard**
- ADJTENS Adjustable Tension Back Upholstery \$142
- L8601 Omit Back Upholstery Less net (\$20)

Arm Heights			
AT99	Range adjusts in 1" increments		8"-10"
CAD			10"-14"
CAF			10"-14"
Range adjusts in 1" increments			
TAD	Low Range	6"-10"	
Range adjusts in 1/2" increments			
TAF	High Range	8"-14"	
Range adjusts in 1" increments			
All measurements from top of seat rail to top of armrest. On AT99 add 1" for upholstered ampads.			

Front Seat-To-Floor Heights					
Fork	Spacer	CASTER			
		AT900 3"	4SPLT 4ALUM2 4"	5SOFT AT901 5ALUM2 5"	AT902 6ALUM2 6"
BRG1 (Short)	None	15.5 16.0	15.5 16.0 16.5	16.5 17.0	17.5
	1"	16.5 17.0	16.5 17.0 17.5	17.5 18.0	18.5
	2"	17.5 18.0	17.5 18.0 18.5	18.5 19.0	19.5
BRG2 (Long)	None	17.5 18.0	17.5 18.0 18.5	18.0 18.5 19.0	18.5 19.0 19.5
	1"	18.5 19.0	18.5 19.0 19.5	19.0 19.5 20.0	19.5 20.0 20.5
	2"	19.5 20.0	19.5 20.0 20.5	20.0 20.5 21.0	20.5 21.0 21.5
SOFR (Suspension)	None	16.5 17.0	16.5 17.0 17.5	17.5 18.0	18.5
	1"	17.5 18.0	17.5 18.0 18.5	18.5 19.0	19.5
	2"	NA	NA	NA	NA

**CASTER SELECTIONS**

**CASTER OPTIONS**

- AT900 3" Roller Blades No Charge  
*Not available with HD - Heavy Duty Option*
- 4SPLT 4" Composite Urethane Light Up No Charge  
*Not available with HD - Heavy Duty Option*
- 4ALUM2 4" Frog Legs® Aluminum Soft Roll Casters (pair) \$134
- 5SOFT 5 x 1 Urethane Casters No Charge
- 5ALUM2 5" Frog Legs Aluminum Soft Roll Casters (pair) \$134
- 6ALUM2 6" Frog Legs Aluminum Soft Roll Casters (pair) \$134
- 41DC 4 in x 1 in Endura Roll Caster (Urethane) \$63  
*Not available with HD - Heavy Duty Option*
- 51DC 5 in x 1 in Endura Roll Caster (Urethane) No Charge  
*Not available with HD - Heavy Duty Option*
- 61DC 6 in x 1 in Endura Roll Caster (Urethane) No Charge  
*Not available with HD - Heavy Duty Option*

**HARDWARE ANODIZING PACKAGE**

*Forks, spacers, caster housing and axle plate are included in the hardware anodizing package. Not available with PROX4F70 or PROX4F80*

- BLKANO Bold Black Anodized Accent Pkg No Charge
- SVRANO Classic Silver Anodized Pkg No Charge

**FORK TYPE (NON SUSPENSION)**

- BRG1 Fork (3 Hole Short Fork) (pair) No Charge
- BRG2 Fork (3 Hole Long Fork) (pair) No Charge

**FORK SPACER**

- SPCR0 No Fork Spacer No Charge
- SPCR1 1" Fork Spacer No Charge
- SPCR2 2" Fork Spacer No Charge  
*Not available with SOFR-F.*

**REAR WHEEL SELECTIONS**

**WHEEL OPTIONS**

- Spoke wheels come standard with dome, quick-release button axle.*
- 20COM 20" Composite No Charge  
*Available with PNEU, PNEUFF or URETH tires.*
- 22COM 22" Composite No Charge  
*Available with PNEU, PNEUFF or URETH tires.*
- 24COM7 24" Composite - 7 Spoke No Charge  
*Available with PNEU, PNEUFF, URETH, or KIKB tires.*
- 559TS 559 Turbo Spoke \$332  
*Not available with HD - Heavy Duty Option. Available with PRIMO or PRIMOTRD.*
- SPIN24 24" Spineray® Spox Wheels \$731  
*Not available with HD - Heavy Duty Option. Available with PRIMO tires.*
- SPIN25 25" Spineray Spox Wheels \$731  
*Not available with HD - Heavy Duty Option. Available with PRIMO or PRIMOTRD.*
- SPIN26 26" Spineray Spox Wheels \$731  
*Not available with HD - Heavy Duty Option. Available with PRIMO tires.*
- 24HPS 24" High Performance Spoke \$210  
*Not available with HD - Heavy Duty Option. Available with PNEU, PNEUFF, PRIMO, or KIKB tires.*
- 26HPS 26" High Performance Spoke \$210  
*Not available with HD - Heavy Duty Option. Available with PNEU or PRIMO tires.*
- SPLX24 24" Spineray® LX Wheels - 12 Spoke \$959  
*Not available with HD - Heavy Duty Option*
- SPLX25 25" Spineray® LX Wheels - 12 Spoke \$959  
*Not available with HD - Heavy Duty Option*
- FN1624 24" Sun Fusion 16 Spoke Wheel \$460  
*Not available with HD - Heavy Duty Option*
- FN1625 25" Sun Fusion 16 Spoke Wheel \$442  
*Not available with HD - Heavy Duty Option*
- FN1626 26" Sun Fusion 16 Spoke Wheel \$442  
*Not available with HD - Heavy Duty Option*
- 20ITESPK 20" Light Weight Perf. Stylized Spoke No Charge  
*Not available with HD - Heavy Duty Option. Available with PNEU, or PNEUFF tires.*
- 22ITESPK 22" Light Weight Perf. Stylized Spoke No Charge  
*Not available with HD - Heavy Duty Option. Available with PNEU, or PNEUFF tires.*
- 24ITESPK 24" Light Weight Perf. Stylized Spoke No Charge  
*Not available with HD - Heavy Duty Option. Available with PNEU, PNEUFF, PRIMO, or KIKB tires.*
- 26ITESPK 26" Light Weight Perf. Stylized Spoke \$105  
*Not available with HD - Heavy Duty Option. Available with PNEU or PRIMO tires.*

*Rear Seat-to-Floor is factory set at 1" - 1.5" lower than Front Seat-to-Floor.*

Rear Seat-To-Floor Ranges	
SEAT-TO-FLOOR HEIGHT	REAR WHEEL
13.5"-18"	20COM, 20SP
14.5"-19"	22COM, 22SP
15.5"-20"	24COM7, 24SP, 24HSP, 24SPIN
16"-20.5"	25SPIN, 559TS, 26SP, 26HPS, 26SPIN

\*Spineray Spoke Color (If color is not chosen, spokes will be black.)

- SPCY - Yellow,  SPCBK - Black,  SPCW - White,  SPCR - Red,  SPCBL - Blue

**TIRE OPTIONS**

<input type="checkbox"/>	PNEU	Pneumatic	No Charge
<input type="checkbox"/>	PNEUFF	Pneumatic - Flat Free	\$128
<input type="checkbox"/>	URETH	Urethane (Treaded)	No Charge
<input type="checkbox"/>	PRIMO	Primo	\$105
<i>Not available with HD - Heavy Duty Option</i>			
<input type="checkbox"/>	PRIMOTRD	Primo 559 Treaded Pneumatic	\$105
<i>Not available with HD - Heavy Duty Option, 25" tire only available with SPIN25, SPLX25 or 559TS.</i>			
<input type="checkbox"/>	KIKB	KIK®-Black	No Charge <b>NEW - No Charge option</b>
<i>Not available with HD - Heavy Duty Option</i>			
<input type="checkbox"/>	CRF123	Schwalbe 1" Blk Flat resistant Non Marking Tire	\$118
<i>The Schwalbe tire is a black high-pressure pneumatic tire that is non-marking, and includes a Smart Guard® puncture resistant layer.</i>			
<i>Schwalbe tire is only available with: 24SP,26SP,24HPS,26HPS,SPIN24,SPIN25,SPIN26.</i>			

**HANDRIM OPTIONS**

<input type="checkbox"/>	ALUMPC	Aluminum Anodized Handrims	<b>Standard</b>
<input type="checkbox"/>	A1512	Black Plastic Coated Handrims	\$114
<input type="checkbox"/>	TIPC	Titanium Handrims	\$651
<i>Not available with HD - Heavy Duty Option. Not available with 20COM, 22COM7, 24COM7, 20SP, 22SP, 559TS.</i>			
<input type="checkbox"/>	NFITSD	Natural Fit Handrims Standard Grip (pair)	\$381
<i>Only available with 24/25/26" spoke wheels.</i>			
<input type="checkbox"/>	NFITSP	Natural Fit Handrims Super Grip (pair)	\$444
<i>Only available with 24/25/26" spoke wheels.</i>			
<input type="checkbox"/>	NFSDLT	Natural Fit LT Handrims Standard Grip	\$412
<i>Only available with 24" composite 7-spoke wheels.</i>			
<input type="checkbox"/>	NFSPLT	Natural Fit LT Handrims Super Grip	\$475
<i>Only available with 24" composite 7-spoke wheels.</i>			

**MOUNTING HANDRIMS (SPOKE WHEELS ONLY)**

<i>Not available with HD - Heavy Duty Option</i>			
<input type="checkbox"/>	IN	Mount Inside Handrim	No Charge
<input type="checkbox"/>	OUT	Mount Outside Handrim	<b>Standard</b>

**AXLES**

<input type="checkbox"/>	A1577	Quick Release Axles	<b>Standard</b>
<input type="checkbox"/>	A1573	Quad Release Axles	\$163
<input type="checkbox"/>	PAXLE	Permanent Axles	No Charge
<i>Standard with Heavy Duty Package.</i>			

**CAMBER ANGLE**

<input type="checkbox"/>	C00	0° Camber	No Charge
<input type="checkbox"/>	C03	3° Camber	No Charge
<input type="checkbox"/>	C06	6° Camber	No Charge

**FRAME FINISH SELECTIONS****FRAME FINISH**

<input type="checkbox"/>	24P	Wet Black	No Charge
<input type="checkbox"/>	29P	White Snow	No Charge
<input type="checkbox"/>	30P	Sunny Yellow	No Charge
<input type="checkbox"/>	60P	Silver Metallic	No Charge
<input type="checkbox"/>	61P	Electric Red	No Charge
<input type="checkbox"/>	62P	Electric Blue	No Charge
<input type="checkbox"/>	105P	Midnight Blue	No Charge
<input type="checkbox"/>	115P	Black Prism	No Charge
<input type="checkbox"/>	116P	Sunkissed Orange	No Charge
<input type="checkbox"/>	117P	Tangerine Red	No Charge
<input type="checkbox"/>	119P	Bubblegum Pink	No Charge
<input type="checkbox"/>	120P	Grape Madness	No Charge
<input type="checkbox"/>	127P	Grasshopper Green	No Charge

**POSITIONING STRAPS/ACCESSORIES****SEAT POSITIONING STRAPS**

<input type="checkbox"/>	A1515	Auto Style Seat Positioning Strap	\$84
<input type="checkbox"/>	A1321	Hook & Loop Seat Positioning Strap	\$66

**BACK ACCESSORIES**

<input type="checkbox"/>	BP02	Back Pack	\$105
--------------------------	------	-----------	-------

**SEAT ACCESSORIES**

<input type="checkbox"/>	SP-01	Seat Pouch	\$83
--------------------------	-------	------------	------

**ARM ACCESSORIES**

<i>Available with AT99 arms only.</i>			
<input type="checkbox"/>	A1571	Fabric Clothing Protectors	\$118
<input type="checkbox"/>	1671-SM	Rigid Side Guards, Small	\$203
<input type="checkbox"/>	1671-LG	Rigid Side Guards, Large	\$203
<input type="checkbox"/>	1671-SMH	Quad Rigid Side Guards, Small	\$203
<input type="checkbox"/>	1671-LGH	Quad Rigid Side Guards, Large	\$203

**PERFORMANCE DETAIL PACKAGE**

- SPORT ProSPIN Sport Performance Detail Package \$70

*Sport Package includes impact guards, rear frame protectors, and clear coat back cane brackets.*

**MISCELLANEOUS ACCESSORIES**

- A1307 Spoke Protectors \$153
- A1580 Rear Anti-Tippers \$109
- 1209 Tire Pump \$63
- A1541 Caster Swivel Locks \$132
- EXTPT Extended Amputee Axle Plate \$114
- Not available with HD - Heavy Duty Option. Only available with 0 degree of camber. Not available with HUMV undermount wheel locks, cannot mount the wheellock far enough back to engage wheels. Not Available with Anti Tippers.*
- 6WLEXT 6" Wheel Lock Extensions \$63

**WHEEL LOCK SELECTIONS**

**WHEEL LOCKS**

- WLPULL Pull To Lock Clamp-On Wheel Locks No Charge
  - WLPUSH Push to Lock Wheel Locks No Charge
  - HUMV Hideaway Undermount Wheel Locks \$149
- Not available with HD - Heavy Duty Option. Not available with URETH tires. Not available with D14, D15, D16 seat depths.*

**FRONT RIGGING SELECTIONS**

**FRONT RIGGINGS**

- 60SSNTAPER 60° Swing-In Swing-Out Footrests (pair) \$241
- 70SSNTAPER 70° Non-Tapered Swing-In Swing-Out Footrests (pair) \$241
- 70SSTAPER 70° Tapered Swing-In Swing-Out Footrests(pair) \$281
- ELRSS Elevating Legrests (pair) \$567

*Not Available with WLPUSH. Interference occurs.*

**FOOTPLATE TYPE (PROX4S)**

- 1651 Flip-Up Composite Footplate No Charge
- Not available with Elevation Legrests.*
- AT5543 Adjustable Angle Flip Up Footplates \$132
- Dorsi/Plantar angle adjustable.*

**FOOTPLATE TYPE (PROX4F70 & PROX4F80)**

- AT5545 One Piece Adjustable Angle Flip-up Footplate \$132
- Not available with HD - Heavy Duty Option*

**FOOTREST/LEGREST ACCESSORIES**

- A1337 Calf Strap \$63
- FRMPRBK ChairArmor Black Neoprene Frame Protectors No Charge
- FRMPRR ChairArmor Red Neoprene Frame Protectors \$70
- FRMPRB ChairArmor Blue Neoprene Frame Protectors \$70
- ALPT Longer Pivot And Slide Tube \$132
- Not available with AT5545 & AT5540.*
- A9342 "H" Heel Strap \$94
- HLOOP Heel Loops without Ankle Straps No Charge

Knee-To-Heel Range (PROX4S)		
Front Rigging Type	Min.	Max.
ELRSS	15.25	21.25
60°	14.00	18.75
70°	13.75	18.75
70° Tapered	14.50	19.50
Range adjusts in 1" increments		
Measurements taken from top of seat rail to back of footrests.		

Knee-To-Heel Range (PROX4F70 & PROX4F80)		
Front Rigging Type	Min.	Max.
Adjustable Angle Flip Up	14.00	18.00
1 Pc. Tubular Flip Up	14.00	18.00
Range adjusts in 1/4" increments		
Measurements taken from top of seat rail to back of footrests.		

WEIGHT LIMITATION 250 POUNDS	
Invacare recommends all Invacare wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.	

**PROSPIN X4 WEIGHT CHARTS**

Frame Weight Swing-Away (lbs)									
Depth	14"	15"	16"	17"	18"	19"	20"	21"	22"
Width									
14"	7.25	7.31	7.37	7.43	7.49				
15"	7.31	7.37	7.44	7.5	7.57				
16"	7.37	7.44	7.5	7.56	7.63				
17"	7.44	7.5	7.56	7.62	7.69				
18"	7.5	7.56	7.62	7.68	7.75	7.81	7.87		
19"			7.68	7.74	7.8	7.87	7.93		
20"			7.74	7.8	7.86	7.93	8	12.42	12.42
22"									12.42

Frame Weight Fixed (lbs)									
Depth	14"	15"	16"	17"	18"	19"	20"	21"	22"
Width									
14"	6.25	6.31	6.37	6.43	6.49				
15"	6.31	6.37	6.44	6.5	6.57				
16"	6.37	6.44	6.5	6.56	6.63				
17"	6.44	6.5	6.56	6.62	6.69				
18"	6.5	6.56	6.62	6.68	6.75	6.81	7.87		
19"			6.68	6.74	6.8	6.87	6.93		
20"			6.74	6.8	6.86	6.93	7	N/A	N/A
22"									N/A



<b>Basic Chair Includes</b>	
•	Folding Curved Frame w/Swingaway Front
•	Seat Width (W14" - W19")
•	Seat Depth (D14"- D19")
•	Swingaway Padded Arms
•	Front Seat To Floor Height Selection (17"-21")
•	Adjustable Height Back Selection w/o Push Handles
•	Slip On Back Upholstery
•	Powder Coated Frame Finish
•	Caster Selection (SPLT4, AT902, AT901)
•	Fork Type (BRG1 or BRG2)
•	Fork Spacer Type (SPCR0, SPCR1, SPCR2)
•	Wheel Options (20SP, 22SP, 24SP, 20COM, 22COM, 24COM7)
•	Urethane (treaded) Tire (URETH)
•	Aluminum Anodized Handrims (ALUMPC)
•	Axle Type (PAXLE or A1577)
•	Wheel Lock Selections (Push or Pull to Lock)
•	Camber Angle Selection
•	Handrim Options (AHDRM or CHDRM)
•	Swingaway Front Rigging
•	Flip Up Composite Footplates

<b>SPECIFICATIONS</b>	
• Frame Style	Folding w/Swingaway Front End
• Seat Widths	14" - 20" & 22"
• Seat Depths	14" - 22"
• Back Height	10" up to 20" (See Back Type for ranges)
• Adjustable Angle Back	0°, 5°, 10°, 15°
• Camber Angles	0°, 3°, 6°, 9°
• Tubing O.D.	1" Aluminum
• Weight Of Chair	15 lbs.* 30 lbs.**
• Weight Limit	250 lbs
Heavy Duty	350 lbs
• Warranty	Lifetime on frame and crossbraces
*16"x16", BACK1216, SLIPBACK, 4ALUM2, BRG1, SPCR1, no rear wheels, no footrests, no anti-tippers, no wheellocks, no arms.	
**16"x16", AT99, BACK1216, SLIPBACK, 4ALUM2, BRG1, SPCR1, SPIN24, PRIMO, ALUMPC, A1577, 70SSNTAPER, 1651, HLOOP, IMPACTG, 1671LG, WLPULL.	

Invacare Canada reserves the right to change prices and specifications without notice.

HEAD OFFICE - 570 Matheson Blvd, Unit# 8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)

Invacare®, and Yes, you can., are trademarks of Invacare Corporation.

[WWW.INVACARE.CA](http://WWW.INVACARE.CA)