



Reset Form

Ultra Low Maxx Power Positioning Accessories
US PRICE LIST AND ORDER FORM

Price Effective August 21, 2019



PLEASE NOTE:

These Accessory Additions are only available with the purchase of a new power positioning system and completed order form TRD0335, TRD0360, TRD0396, TRD0489, TRD0498, TRD0506 and TRD0507 with all Client Specifications completed. Subject to approval

Dealer Name:
Dealer Account Number:
PO #:
Tag:

LATERAL TILT MODULE

Power Lateral Tilt Module(up to 15°) Part # LATM... \$3,995
All options must be selected on the system order form prior to selecting this option.
Reduced Drive while in Lateral Tilt? Yes Part # LATDY No Part # LATDN

NOTES:

- 1. Weight capacity is 250 lbs.
2. Only available on TDX-SP2, 3G Storm, TDX-SP, 6000ZHD, Edge HD and Rovi A3 Bases.
3. When ordering a lateral tilt module on a Storm short base frame, standard fixed caster must be ordered on the base.
4. Only available on 16" to 20" seat widths and 16" to 22" seat depths.
5. Only available with Multi-function electronics Part # MFCBD.
6. Must order the Group 24 battery tray for the TDX-SP2 and TDX-SP, which will utilize Group 22 batteries to accommodate the lateral tilt module.

POWER ANTERIOR TILT MODULE

Power Anterior Tilt Module(up to 45°) Part # PAT45 . \$3,995
All options must be selected on the system order form prior to selecting this option.
Synchronized Option for Center Mount Front Rigging with Anterior Tilt Part # SOCM \$1,350
For LNX, Maxx Style Fixed and Power Elevating Calf Panel.
Synchronized Option Recline with Anterior Tilt Part # SOR \$500
For use on Tilt/Recline Systems

NOTES:

- 1. This module may limit front rigging options.
2. Weight capacity 250 lbs.
3. Only available on 16" to 20" seat widths and 16" to 22" seat depths.
4. Only available with Multi-Function electronics Part # MFCBD.
5. All front rigging will be set a 70° unless using the Synchronized Center Mount option.
6. Minimum 18" seat to floor height at zero degrees.
7. Minimum finished knee to heel to touch floor is 13".
8. This is an upgradeable option only.

UPFRONT MODULE **ROVI BASES ONLY**

Upfront Module (up to 45°) with Anterior Tilt and 7" of Seat Elevation Part # UM... \$5,995
(System moves Forward 5" while elevating forward and includes a synchronized leg mechanism to tuck when elevated forward and when extends forward in anterior tilt)
All options must be selected on the system order form prior to selecting this option.
Reduced Drive while Upfront? Yes Part # UMDY No Part # UMDN

NOTES:

- 1. Weight capacity is 250 lbs.
2. Only available with LNX Power Center Mount Front Rigging or Maxx Style Manual Center Mount foot platform.
3. Recline is recommended to assist with forward transfer.
4. Only available on 16" to 20" seat widths and 16" to 22" seat depths.
5. Only available with Multi-Function electronics Part # MFCBD.
6. Articulating Vents and O2 Holders are not available.
7. Power Sliding Back (PSB) is not available on the upfront module.
8. Power Telescoping and Power Flip up foot platforms are not available.

POWER SWING AWAY JOYSTICKS

Power Swing Away Joystick Part # JPO... \$2,995
(Joystick Power "on" operation)
Select Movement: Horizontal (Quad link style) Part # PSAJH Vertical Flip up (Mid line style) Part # PSAJV
Select Mounting Position: Left Side Right Side Part # ML MR

Operate Power Joystick/Chin Boom with Joystick Power "off" option Part # JPO... \$750

NOTES:

- 1. Comes with an egg switch
2. Counts as an additional electronic function.
3. Allows power joystick to operate when chair is turned off.

POWER SWING AWAY CHIN BOOM

Power Swing Chin Boom Part # PSACB... \$2,995
(Joystick Power "on" operation)
Select Mounting Position: Left Side Right Side Part # ML MR

Operate Power Joystick/Chin Boom with Joystick Power "off" option Part # JPO... \$750

NOTES:

- 1. Comes with an egg switch
2. Counts as an additional electronic function.
3. Allows power chin boom to operate when chair is turned off.

TRANSFER ARM

Power Transfer Armrest Part # PTAL... \$1,995
Select Mounting Position: Left Side Right Side Part # PTAL PTAR

Manual Transfer Armrest Part # MTAL... \$995
Select Mounting Position: Left Side Right Side Part # MTAL MTAR

Power Flip Back Armrest Part # PFBAL... \$1,995
Select Mounting Position: Left Side Right Side Part # PFBAL PFBAR

NOTES:

- 1. Weight capacity is 200 lbs.
2. Only with Adjust Dual Post Style Armrest for Tilt/Recline, arm height 9"-13"
3. Required Multi-function electronics ordered on chair Part # MFCBD.
4. Minimum seat depth available is 16" and 16"-20" seat width.
5. Minimum seat to floor height 18"
6. Transfer only available in the chair home position.
7. Not available with Upfront Module.

Ultra Low Maxx Power Positioning Accessories

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

POWER COMMUNICATION DEVICE (1)(2)(3)(4)(5)(6)

Power Communication Device - Single Actuator⁽¹⁾⁽²⁾ \$2,995
(Joystick Power "on" operation)

Select Mounting Position: Left Side Right Side
Part # PCDSL PCDSR

Power Communication Device - Dual Actuator⁽¹⁾⁽²⁾ \$3,995
(Joystick Power "on" operation)

Select Mounting Position: Left Side Right Side
Part # PCDDL PCDDR

Operate Power Joystick/Chin Boom with Joystick Power "off" option⁽³⁾ Part # JPO \$750

Additional Electronics⁽³⁾ Part # AE \$3,995

NOTES:

1. Considered an additional electronic function and requires Multi-function electronics.
2. Consider as an additional two functions.
3. May require additional electronics if system has more than 6 actuators (functions).
4. Minimum seat depth available is 16".
5. Weight capacity 20 lbs.
6. Comes with an egg switch.
7. Allows Power Communication Device to operate when chair is turned off.

POWER FOOT PLATFORM OPTIONS FOR LNX

For use with LNX Foot platform ordered separately

Power Telescoping Foot Platform for LNX⁽¹⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾
Powers lower leg length vertically up to a maximum of 4" extension.

Part # LNXPTF **\$2,250**

Power Flip up Foot Platform for LNX⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾⁽⁶⁾
Part # LNXPFUF **\$1,995**

NOTES:

1. Available with one piece calf panel only and not available with a finished knee-to-heel of 13" or less. Requires 2-6 function electronics. Not available with Split Footplates.
2. Comes standard with a 12" x 10.75" footplate made of durable aluminum 1/4" thick with heel loops.
3. The front rigging will be set at approximately 80° on the Storm 3G.
4. Minimum finished knee to heel is 13".
5. Minimum seat to floor height is 17".
6. Not available on Upfront Module.

POWER FOOT PLATFORM OPTIONS

Fixed Center Mount footplatform is included

Power Telescoping Foot Platform with Fixed Center Mount⁽¹⁾⁽³⁾⁽⁴⁾⁽⁵⁾ Part # FCMFPTF \$2,845
Powers lower leg length vertically up to a maximum of 4" extension.

Power Flip up Foot Platform with Fixed Center Mount⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾⁽⁵⁾ Part # FCMFPFUF . . . \$2,590

NOTES:

1. Included the Fixed Center Mount Foot platform which does not elevate.
2. Comes standard with a 12" x 10.75" footplate made of durable aluminum 1/4" thick with heel loops. Optional sizes available per standard order form.
3. This option is not available Storm 3G unless Flared Front Casters for Storm short base are added.
4. Minimum finished knee to heel is 13".
5. Minimum seat to floor height is 17".

POWER CALF PANEL (1)(2)(3)(4)


Power Elevating Calf Panel Part # PECP \$3,995

NOTES:

1. No footplate available.
2. Calf panel measures up to 15" in length and measures 13" wide at the bottom and tapers up to within 2" narrower than the seat width.
3. Minimum seat depth 16" deep.
4. Not available for use with Dual Front Rigging adaptor.

CUSTOMIZED FOOT PLATFORM OPTIONS

Custom Foot Platform⁽¹⁾ Part # CFP \$550



Onyx Pad for Any Foot Platform⁽²⁾⁽³⁾ Part # OCUFP. . . . \$450

NOTES:

1. Available for Maxx/LNX Center Mount.
2. Available for all Maxx/LNX Center Mount Foot Platform options.
3. Pad thickness is 1"

CUSTOMIZED FOOT PLATFORM OPTIONS⁽¹⁾⁽²⁾⁽³⁾⁽⁴⁾

Power Assisted Legrest Swing Away Mechanism . \$3,995

Specify operated through Switch Joystick
Part # PALSAMS PALSAMJ

NOTES:

1. Only available with Multi-Function electronics
2. Compatible with Motion HD 70° Swing Away footrest, Maxx Style Swing Away Legrest Power or Manual.
3. Legrests are not removable.
4. Minimum 16" seat width and minimum 16" seat depth.

STORM FLARED CASTER OPTION

Flared Front Casters for Storm Short Base⁽¹⁾

Part # FFC **\$950**

NOTE:

1. Allows 90° front rigging on Short Storm bases, comes with 10" x 10" footplate.

OTHER ACCESSORIES⁽¹⁾

Power Back Pack Holder with Single Back Pack . \$2,995

Select Mounting Position: Left Side Right Side
Part # BPHL BPHR

Concepts Back Pack Single Part # SBP \$125

Double Part # DBP \$195

NOTE:

1. Weight limit 20lbs.

Ultra Low Maxx Power Positioning Accessories

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

SPECIALTY PAD OPTIONS⁽¹⁾⁽²⁾

Optional Sized Flat One Piece Calf Panel for LNX
 Part # OSCPL **\$450**
 Maximum pad size is seat width less 2" and 10" in length. With minimum 13" finished knee to heel.

Specify preferred size each Height: _____ Width: _____

Lateral Leg Supports for LNX One Piece Calf Panel. . \$650
Select Pad Size::
 XS(5.5"W x 4"H) Part # LLSLNXX S(7"W x 6"H) Part # LLSLNXS
 M(9"W x 4"H) Part # LLSLNXM L(13"W x 4"H) Part # LLSLNLX

Lateral Leg Supports for Power Elevating Calf Panel \$650
Select Pad Size::
 XS(5.5"W x 4"H) Part # LLSPCPX S(7"W x 6"H) Part # LLSPCPS
 M(9"W x 4"H) Part # LLSPCPM L(13"W x 4"H) Part # LLSPCPL

Optional Sized Individual Calf Pads for Front Rigging⁽³⁾
 Part # OSICP **\$650**

Specify preferred size each Height: _____ Width: _____

Optional Sized Flat Arm Pads
 Part # OSFAP **\$650**

Specify preferred size each Length: _____ Width: _____

Optional Sized Lateral Pad or Elbow Block
 Part # OSLPEB **\$650**

Specify preferred size each Length: _____ Width: _____

- NOTES:**
1. Subject to clearance available.
 2. Startex fabric over polyurethane foam is standard.
 3. Available on Motion Concepts front rigging options.

ADDITIONAL BATTERY TRAYS

Storm Long Base Additional Front Mounted Battery Tray Group 40 Battery Part # SLBT40 **\$500**

Storm Short Base Additional Front Mounted Battery Tray for U1 Battery Part # SLBTU1 **\$500**

Storm Base 3 x Group 22 Battery Tray Conversion⁽³⁾
 Part # S322 **\$950**

TDX-SP Base Rear Caster Mounted U1 Battery Tray⁽²⁾
 Part # TDXU1 **\$500**

Edge HD 3 x Group 22 Battery Tray Conversion⁽¹⁾
 Part # EDGE322 **\$500**

- NOTES:**
1. Comes with steel battery door cover. Must be ordered with . Group 24 Battery Tray - Wide Base TDX-SP.
 2. Available on Wide Base TDX-SP Only and must have factory . Transport brackets - TBKTS.
 3. Deletes factory battery boxes.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

ADDITIONAL ELECTRONIC OPTIONS⁽¹⁾

Step Down Converter⁽¹⁾ Part # SDC **\$500**
Converts 24v to 12v.

Richardson Powertech Vent Power Center⁽¹⁾
 Part # RPVPC **\$1,950**

- NOTE:**
1. Should this product be used in life support applications Motion Concepts will not assume any risk and shall be held . harmless against all damages.

MEMORY SEATING OPTIONS⁽¹⁾⁽²⁾⁽³⁾

Single Function Memory Seating Part # SFMS **\$495**
**Please call for specification sheet 1-888-433-6818*

Multiple Function Memory Seating Part # MFMS **\$995**
**Please call for specification sheet 1-888-433-6818*

- NOTES:**
1. Does not operate elevation function, anterior tilt function, lateral tilt function or upfront module.
 2. Comes with a quad button switch with LED indicator.
 3. Only available with Multi-function electronics.

SAFETY LIGHTS

Front Safety Lights (pair)
 Operated through switches⁽¹⁾ Part # FSLS **\$495**
 Operated through Joystick⁽²⁾ Part # FSLJ **\$995**

Rear Safety Lights (pair)
 Operated through switches⁽¹⁾ Part # RSLS **\$495**
 Operated through Joystick⁽²⁾ Part # RSLJ **\$995**

Front & Rear Light Combo Package Credit*
 Part # LPC **(\$150)**
**Purchase both sets and receive a credit*

Ground Effects Lighting through Switch
 Part # GEL **\$350**

- NOTES:**
1. Operated through a separate power supply and single push . button switch.
 2. Only available with Multi-Function electronics

Total Retail Value \$ _____

Email Form