



# Freedom<sup>®</sup> NXT<sup>™</sup>

Generation NEXT RIGID Tilt-in-Space Mobility

## Freedom: for Kids (and other people too)

NXT RIGID: Almost twins...

Now, everything you love about the NXT in a rigid frame design!



**RIGID:**  
Lightweight  
& dependable

**RIGID:**  
Strong  
& durable

**RIGID:**  
Easy to grow  
& to use



**Looks can be deceiving...**



**NXT Folding**

Common with twins, the NXT RIGID (left) is almost identical to the NXT folder (above). It has all the same great features and options.

Both are available in three pediatric and three adult sizes. And, at 23 lbs., the NXT RIGID is the lightest, easiest to use, and most durable tilt-in-space available. You'll be happy you roll with FREEDOM.

**Innovative "firsts" to keep kids safe**

The NXT is the first tilt-in-space wheelchair to have an attendant:

- drum braking system to help control speed, direction, and lock the wheels
- hand control that both brakes and locks the wheels and tilts the frame
- anti-tip with a built-in curb assist
- plus, the NXT incorporates seat frame leveling guides so you know the frame is straight and true when not in tilt.



**NXT RIGID**



**NXT Future Fit**

When the growth spurts hit and the frame won't fit, it's Freedom to the rescue. Freedom's "NXT Future Fit" program is a one time, no cost program to provide parts to make the frame wider. Freedom's NXT is a good fit today, tomorrow and forever!

**The NXT is versatile, quickly transforming for ever-changing kids**

**Freedom® NXT™: Generation NEXT RIGID Tilt-in-Space Mobility Specifications**

**Seat Frame Widths**

The NXT is available in 11 sizes from 10" wide to 20" wide.  
 Each size from 10" wide to 16" wide adjusts an additional 2" in width.  
 The 17" wide size adjusts an additional 1" in width to 18" wide.  
 The 18", 19" and 20" wide sizes are set.

**Frame Color**

Available in **25** colors.

See NXT order form for further details.

**Seat Frame Depths**

The NXT is available in 3 sizes: short, medium or long.  
 Each size adjusts an additional 4" in depth.  
 Seat frame depth is measured without seating, from the front of the back post to the front end of the seat frame.  
 The short range is 12" to 16."  
 The medium range is 14" to 18."  
 The long range is 16" to 20."

**Seat Depth Extensions**

For additional seat frame depths, extensions can be added to the back of the seat frame. There is a 2", a 3" or a 4" extension to choose from.

**Transport Package**

Includes four transportation tie down securement points attached to the NXT frame. Package must be specified "on chair" at time of order. Available for pediatric and adult sized NXT mobility bases.

**Seat to Floor Height**

The NXT can be configured in a variety of seat to floor heights.  
 Selections in front casters and rear wheel sizes determine a range of choices from 12.5" low to 18" high, and most everything in between.

**Tilt Range**

45 degrees.

**Rider Capacity**

225 lbs.



*Our kids are our passion.™*

**Freedom Designs Incorporated**  
 2241 N. Madera Road, Simi Valley, CA 93065

phone: 800-331-8551  
 fax: 888-582-1509  
 e: CustomerService@freedomdesigns.com  
 www.freedomdesigns.com

Freedom Designs Incorporated (FDI) maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without notification. Check with FDI for the latest information.



Copyright © 2014 Freedom Designs, Inc. All rights reserved. Printed in the USA. WBNXTR 112014



**Manufactured in the U.S.A.**