



Reset Form



REQUIRED INFORMATION

Date of Order: _____ Dealer Account #: _____
 Dealer Name: _____
 PO #: _____ Tag _____
 After Market Items shipping to Motion Concepts: Yes No
 List Items: _____
 PO #: _____
 Name of the Manufacturer: _____
 Purchasing Contact: _____
 Phone: _____ Fax: _____
 E-mail: _____
 RTS/Therapist: _____
 Ship to Address: _____
 City: _____ State: _____
 Zip: _____ Client Gender: M F
 Special Client Conditions: _____

REQUIRED MEASUREMENTS

To ensure system is accurately configured please fill in all required measurements.

Client Height:	
Client Weight (lbs):	
A. Seat to Shoulder:	
B. Trunk Depth:	
C. Chest Width:	
D. Knee to Back:	
E. Seat to Top of Head:	
F. Elbow to Hand:	
G. Seat to Elbow:	
H. Hip Width:	
I. Knee to Heel:	
Cushion Thickness:	

Weight Capacity is up to 300lbs unless specified otherwise



HCPCS CODES	K0848, K0849
SUSPENSION SYSTEM	Active Ride Control (ARC)
BASE WIDTH	25.5"
BASE LENGTH	36.5"
DRIVE WHEEL SIZE ⁽¹⁾	14"
CASTER WHEEL SIZE (front/rear)	6"
BASE ONLY TURNING RADIUS	20.5"
DRIVETRAIN	4-pole EAD Motors

BATTERIES	2 x M34 3 x M50 option for Vents
BATTERY CHARGER	8A, Off Board
BATTERIES WEIGHT (each)	38.5 lbs
RANGE ⁽²⁾	17.7 miles
WEIGHT CAPACITY	300 lbs
BASE WEIGHT (w/o batteries)	157 lbs
GROUND CLEARANCE	3"
MAXIMUM SPEED	6.2 mph

BASE VITALS

⁽¹⁾ All A3 wheels come with solid inserts
⁽²⁾ Actual driving range and speed may vary due to factors such as user weight, type and grade of terrain, battery condition & charge level; type and condition of the drive wheels, and variations in the drive and control systems.

SOLID SEAT PAN - K0848

CAPTAIN SEAT - K0849



Contoured Static Seating on ROVI A3 Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

BASE COLOR

Select Color: Standard With Lights⁽¹⁾ \$1,195

- | | | |
|--|---|---|
| <input type="checkbox"/> Blizzard White | <input type="checkbox"/> Part # 800914-00200-03 | <input type="checkbox"/> Part # 800914-00201-03 |
| <input type="checkbox"/> Matte Black | <input type="checkbox"/> Part # 800914-00200-50 | <input type="checkbox"/> Part # 800914-00201-50 |
| <input type="checkbox"/> Solar Yellow | <input type="checkbox"/> Part # 800914-00200-69 | <input type="checkbox"/> Part # 800914-00201-69 |
| <input type="checkbox"/> Motion Green | <input type="checkbox"/> Part # 800914-00200-36 | <input type="checkbox"/> Part # 800914-00201-36 |
| <input type="checkbox"/> Original Orange | <input type="checkbox"/> Part # 800914-00200-23 | <input type="checkbox"/> Part # 800914-00201-23 |
| <input type="checkbox"/> Bright Blue | <input type="checkbox"/> Part # 800914-00200-40 | <input type="checkbox"/> Part # 800914-00201-40 |
| <input type="checkbox"/> Red Delicious | <input type="checkbox"/> Part # 800914-00200-21 | <input type="checkbox"/> Part # 800914-00201-21 |
| <input type="checkbox"/> Party Pink | <input type="checkbox"/> Part # 800914-00200-2C | <input type="checkbox"/> Part # 800914-00201-2C |
| <input type="checkbox"/> Purple Reign | <input type="checkbox"/> Part # 800914-00200-80 | <input type="checkbox"/> Part # 800914-00201-80 |
- Matching Battery Doors (select fender color first) . . \$495

NOTES:

1. Fully integrated LED lighting kit. Includes headlights, taillights and turn signals. Must select Joystick for Lights feature or Omni for operation.



Please see Custom Design Options on page 5/6

TIRE OPTIONS

- | | |
|---|--------------|
| <input type="checkbox"/> Standard Solid Tires Part # STIRES | STD |
| <input type="checkbox"/> Pneumatic Tires Part # PTIRES | \$199 |

MILITARY DECAL

3" Diameter Decal (Uninstalled) \$50



- | | |
|------------------------------------|----------------------------------|
| <input type="checkbox"/> Army | <input type="checkbox"/> Navy |
| <input type="checkbox"/> Air Force | <input type="checkbox"/> Marines |



Please note that the HCPCS codes indicated on this order form are our interpretation of the new codes and how we believe our products should be coded. The individual DMERCs will determine if our interpretation is appropriate and whether or not items coded will be accepted or denied. Motion Concepts cannot guarantee that use of the designated codes will ensure coverage or payment. Contact the insurer to confirm compliance with coding and coverage policy

FULL CAPTAIN-STYLE SEAT - K0849

All bases below come standard with Full Captain-Style Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

- ROVI A3 5.8 mph Quick Speed Motors** ⁽¹⁾
 HCPCS Code: K0849 Part # 888WNLEA3-MP6-CS **\$7,095**
- Optional 6.2 mph High Speed Motors**⁽²⁾ **\$950**
 with OnTraxx Enhanced Driving Module Part # HSM62
- OnTraxx Enhanced Driving Module**
 HCPCS Code: K0108 Part # OEDM. **\$1,550**

- Seat Size: Grey** 18"x18" Part # SV18171/222RG
Black 18"x18" Part # SV18171/222RB

STF for 18 x 18⁽²⁾: 19" 20" 21" 22"

- Seat Size: Grey** 20"x20" Part # SV2121211/2RG
Black 20"x20" Part # SV2121211/2RB

STF for 20 x 20⁽²⁾: 20" 21" 22" 23"

NOTES:

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt, Weight capacity is 300 lbs.
2. Seat to floor height is measured at the seam of cushion on 20" widths and at the top of the leg trough on 18" widths. Seat to floor heights may vary +/- .25" dependent on base conditions.

SOLID SEAT PAN & CAPTAIN-STYLE SEAT - K0848

All bases above come standard with a Solid Seat Pan style Captain Seat in Grey or Black Vinyl and 22" high backrest. Includes a Semi-Recline hinge of 120° with left side only adjustable lever. Standard adjustable height, flip-back armrests and a foot platform that measures 12" wide by 9" deep.

- ROVI A3 5.8 mph Quick Speed Motors** ⁽¹⁾
 HCPCS Code: K0848 Part # 888WNLEA3-MP6-SS **\$6,950**
- Optional 6.2 mph High Speed Motors**⁽²⁾ **\$950**
 with OnTraxx Enhanced Driving Module Part # HSM62
- OnTraxx Enhanced Driving Module**
 HCPCS Code: K0108 Part # OEDM. **\$1,550**

- Select:** **Grey** Seat Size: 18"x18" Part # SP181821RG
 20"x20" Part # SP202023RG
 Black Seat Size: 18"x18" Part # SP181821RB
 20"x20" Part # SP202023RB

Select STF⁽²⁾: 17" 18" 19" 20"

NOTES:

1. Price of the powerbase includes: Active Ride Control suspension, tie-down hardware for 4-point strap system, 8 amp off-board charger, flat free inserts and positioning belt, Weight capacity is 300 lbs.
2. Seat to floor height is measured at the front of the seat pan based on a standard seat size. Seat to floor heights may vary +/- .25" dependent on conditions of the power base and on seat depths greater than 18".

Contoured Static Seating on ROVI A3 Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.


Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

SEAT CUSHION REFERENCE GUIDE



Cushion	Meshtex	Startex	Infection	Spacetex	Fusion	Onyx
Libra	N/A	STD	\$50	N/C	N/A	N/A
PS	N/A	STD - Reversible	STD - Reversible	N/A	N/A	N/A
PSP	N/A	STD - Reversible	STD - Reversible	N/A	N/A	N/A
Vi	STD	\$50	\$50	N/A	\$150	\$100

SEAT CUSHION OPTIONS - FOR SOLID SEAT

- Matrx Libra Seat Cushion 18" w-20" w** HCPCS Code E2624 **\$.485**
- Matrx PS Seat Cushion 18" w - 20" w** HCPCS Code E2605 **\$.345**
- Matrx PSP Seat Cushion 18" w - 20" w** HCPCS Code E2607 **\$.395**
- Matrx Vi Seat Cushion 18" w - 20" w** HCPCS Code E2607 **\$.425**

FABRIC OPTIONS FOR CUSHIONS

- Spacetex Fabric Option for Libra Cushion Part # SFLF N/C
- Infection Control Fabric Option for Libra Cushion
Part # ICFL \$50
- Fusion Fabric Upgrade for Matrx Vi Cushion Part # FFV . . . \$150
- Onyx Fabric Upgrade for Matrx Vi Cushion Part # OCUMA . . . \$100
- Startex Fabric Upgrade for Matrx Vi Cushion
- Smooth Side up Part # SFVS \$50
- Fabric Side up Part # SFVF \$50

SINGLE FUNCTION CONTROL OPTIONS

- R-Net LED Remote Joystick¹** Part # D51122SC \$295
- R-Net LED Remote Joystick for Lights⁽¹⁾**
Part # D51108 \$500
- R-Net Color Remote Joystick with Monojacks⁽¹⁾⁽²⁾**
Part # D51165 \$600
- R-Net Color Remote Joystick with Monojacks for Lights⁽¹⁾** Part # D50982 \$650
- R-Net Advanced Joystick⁽¹⁾** Part # D51711. \$1295
3.5" high-resolution VGA color display. IR Control as standard
Includes Integrated Bluetooth module.
- R-Net Expandable Controller⁽³⁾** \$615
HCPCS Code: E2377 Part # D51109
- R-Net Harness for Expandable Controller⁽³⁾** \$415
HCPCS Code: E2313 Part # REH

- NOTES:**
1. Requires an R-Net Controller and Harness
 2. Light package not available with the joystick.
 3. Required for R-Net Remote Joysticks.

JOYSTICK MOUNT

Fixed for Solid Pan/Captian Seat Part # 1343-00707-00 . . . N/C

Select Mounting Position: Left Right
Part # ML Part # MR

Swing Away Mount Bracket \$295

Select Mounting Position: Left Part # 1343-00709-00-A3
 Right Part # SAQR

JOYSTICK HANDLES

- "T" Handle for Joystick Knob** \$115
HCPCS Code: E2323 Part # SA80941
- Ball Handle for Joystick Knob** \$115
HCPCS Code: E2323 Part # SA72441
- BodyPoint Joystick Handle** HCPCS Code: E2323 \$120
 - 3" U-Shaped Part # PC101B-20
 - 4" U-Shaped Part # PC102B-20
 - Domed Rubber Part # PC107B-20

ELECTRONIC OPTIONS

- R-Net Output Module(EMC)⁽¹⁾⁽²⁾** \$995
Part # D50883
- R-Net Bluetooth Module⁽¹⁾⁽²⁾** \$895
Select OS: iOS Apple Part # D51580 Android/PC Part # D51111
- Smart Phone Holder** Part # SPH N/C

- NOTES:**
1. Requires an R-Net Controller and Harness order above
 2. Must order R-Net Color Remote joystick with these options

ARMRESTS (MUST BE SELECTED)

- Adjustable Height Flip Back Arm Left** Part # 104300-00112. . \$125
- Adjustable Height Flip Back Arm Right** Part # 104400-00112 \$125

OTHER OPTIONS

- Wheelchair Transport Brackets⁽¹⁾** Part # XBKTS \$250

- NOTES:**
1. For unoccupied use only

BATTERIES

- 2 x MK GROUP 34 Batteries** Part # XBKTS M34 SLD G(Two) \$650
HCPCS Code: E2359

Total Retail Value \$ _____

Email Form

Contoured Static Seating on ROVI A3 Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CUSTOM FENDER DESIGN BY BOB RICE

Custom Fender Design by Bob Rice \$500

Standard Fender Configuration Part # Below N/C

Fenders w/Lights Configuration Part # Below + "L" . . . \$1195



Flames
Part # CFF



Personalized
Part # CFP



Camouflage
Part # CFC



Silver Concrete
Part # CFSC



U.S. Navy
Part # CFUSN



U.S. Coast Guard
Part # CFUSCG



U.S. Air Force
Part # CFUSAF



U.S. Marines
Part # CFUSAM



U.S. Army
Part # CFUSAA

DESIGN DESCRIPTION

(PLEASE USE IF YOU SELECTED PERSONALIZED OPTION)

Empty box for design description.

DEALER SIGN OFF

This concept is agreed to and accepted by:

Dealer Name: _____

PO #: _____

Approved with no changes

Approved with changes

Review again

Date: _____

Signature: _____

Contoured Static Seating on ROVI A3 Base

Specific Client Measurements Required - Please call customer service at 1-888-433-6818 and/or fax 1-888-433-6834 to ensure that we have all the client details we require for this type of option.

Dealer Name: _____

Dealer Account Number: _____

PO #: _____

Tag: _____

CUSTOM DIPPED DESIGN OPTIONS

Please select Fender option first, if battery doors are to match the fenders. If not make the design selection below.

Dipped Fender Designs



Standard Fender Configuration
Part # CDF + Code \$395

Fenders w/Lights Configuration
Part # CDFL + Code \$1,590

Front Battery Door Part # FBD + Code \$395

Rear Battery Door Part # RBD + Code \$395

Front/Rear Battery Doors
Combo Credit Part # FRBDC (\$40)

Fenders & Front/Rear Battery Doors
Combo Credit Part # FFRBDC (\$85)



CUSTOM DIPPED DESIGN OPTIONS



U.S. Flag
Code # 1



Camouflage
Code # 2



Electric Circuit
Code # 3



Flames
Code # 4



80's Green Rock
Code # 5



80's Black Rock
Code # 6



Pink Camo
Code # 7



Flower Power
Code # 8



Graffiti
Code # 9



Animal Paw Print
Code # 10