# **INVACARE®**

# MyOn™ SERIES ULTRALIGHTWEIGHT WHEELCHAIRS





The **Invacare® MyOn™ Series Wheelchairs** are the newest ultralightweight folding wheelchairs in the Invacare product family. It matches modern styling and versatility with overall exceptional value.

The Invacare® MyOn™ Wheelchair offers the flexibility to meet a wide variety of needs and preferences. This highly adjustable chair, from the height of the backrest to the armrests, allows for maximum customization and adaptability. Chair adjustments are made with just a few simple tools: 2 Allen keys and a wrench. A tool kit is provided with each MyOn Wheelchair when ordered. Reducing time and complexity is an important factor when adjusting a wheelchair, and the MyOn Wheelchair was designed with this in mind. This chair can easily be adapted for a precise, individualized fit in virtually no time at all.

#### **Features**

Crossbraces for easy folding
Transfer-friendly wheel locks
Fold down back canes
Tire options of Kenda®, Treaded pneumatic
and pneumatic flat free

#### **HCPCS Code: K0005**










- 1. Compact folded footprint
- 2. Standard wheel locks fold down easily for side transfer
- 3. Fold down back canes
- 4. Crossbraces with lack of round tubing for easy folding





<sup>\*16&</sup>quot; x16", 16" Back Height, 5" Composite Hub Casters, 24" Rear Wheels w/ Treaded Pneumatic Tires & Aluminum Anodized Handrims, Padded Swing-Away Arms, Plastic Removable Side Guards, Push-to-Lock Wheel Locks, No Seat Cushion and Anti-tip assembly.

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

The Invacare® MyOn™ Jr. Wheelchair is an excellent choice for children who need a chair with essential adjustments to meet their specific needs. The MyOn Jr. Wheelchair has been developed to match the individual needs of a child and can grow as they grow. Because a child's needs will change over time, the MyOn Jr. Wheelchair offers a large range of options to accompany the child in their development and to provide the necessary clinical support. The MyOn Jr. Wheelchair comes in attractive, modern styles, designed to appeal to kids.

#### **Features**

Kid-friendly spoke protectors Bear shaped side guards Telescoping crossbraces Width and depth adjustable upholstery Tire options of Kenda®, Treaded pneumatic and Pneumatic flat free

#### **HCPCS Code: E1238**













- 1. Compact folded footprint
- 2. Depth adjustable upholstery
- 3. Standard wheel locks fold down easily for side transfer
- 4. Fold down back canes
- 5. Kid-friendly spoke protectors
- 6. Pediatric widths and depths











Back



Front Frame





Model # **FXMYONJR GRMYONJR** 



Seat Dimensions 8"- 15"W 8"- 15"D











Diameter













Knee-to-Heel Range 6.75"- 15"

Center-of-Gravity -1" to +2"

Wheels 20", 22", 24"

Casters 5", 6"

Total Weight\* 27 lb.

Weight Capacity 175 lb.

Limited Warranty 5 Years Frame/

Crossbraces

<sup>\*</sup>Includes: 15"x15", 16" back height, 5" casters, 24" rear wheels with urethane tires and aluminum handrims, desk length T-arms, adjustable width and height nylon back upholstery and push-to-lock wheel locks.

## **Invacare Corporation**

www.invacare.com

#### USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

#### Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

### ©2012

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols <sup>™</sup> and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Kenda is a registered trademark of Kenda Rubber Industrial Co., Ltd. Form No. 12-502 Rev. 12/12 120604

