



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Invacare® Formula® CG with MK6™ for TDX® Spree/SC Retrofit Kit

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWER SEATING BASE

- TDXSPRESEAT Seat for Formula CG System \$755.00
 TDXSCSEAT Seat for Formula CG System \$755.00

SYSTEM TYPE

- CGTK Tilt Only to existing ADJASBA Seat \$5,355.00

SEAT SELECTION

- CTB Contoura™ Seating \$625.00 E2620 E2621
Available on TDX SC only.
 OMCU Omit Contoura Seat Cushion ... Less \$20.00
 CTBMSH Contoura™ seating with mesh upholstery \$635.00
 OMCU Omit Contoura Seat Cushion ... Less \$20.00 E2620 E2621
 BPO Back Posts Only (Round Canes) No Charge
 CURV Curved No Charge
Available on TDX SC only.
 CNB Conventional No Charge
Available on TDX Spree only.

ADJUSTABLE SEAT WIDTH RANGE

- Clear
 1216W 12" to 16" Seat Width Range No Charge
Available on TDX Spree only.
 1620W 16" to 20" Seat Width Range No Charge
Available on TDX SC only.

SEAT WIDTH SETTING

- Clear
 W12 12" Seat Width No Charge
 W13 13" Seat Width No Charge
 W14 14" Seat Width No Charge
 W15 15" Seat Depth No Charge
 W16 16" Seat Width No Charge
 W17 17" Seat Width No Charge
 W18 18" Seat Width No Charge
 W19 19" Seat Width No Charge
 W20 20" Seat Width No Charge

ADJUSTABLE SEAT DEPTH RANGE

- Clear
 1218D 12" to 18" Seat Depth Range No Charge
 1619D 16" to 19" Seat Depth Range No Charge
 1922D 19" to 22" Seat Depth Range No Charge

CONTOURA

- Clear
 1618D 12" to 18" Adjustable No Charge
 1920D 19" to 20" Adjustable No Charge

SEAT DEPTH SETTING

- Clear
 D12 12" Seat Depth No Charge
 D13 13" Seat Depth No Charge
 D14 14" Seat Depth No Charge
 D15 15" Seat Depth No Charge
 D16 16" Seat Depth No Charge
 D17 17" Seat Depth No Charge
 D18 18" Seat Depth No Charge
 D19 19" Seat Depth No Charge
 D20 20" Seat Depth No Charge

BACK HEIGHT

BACK HEIGHT

- Clear
 TSBH16 16" Back Height No Charge
 TSBH17 17" Back Height No Charge
 TSBH18 18" Back Height No Charge
 TSBH19 19" Back Height No Charge
 TSBH20 20" Back Height No Charge
 TSBH21 21" Back Height No Charge
 TSBH22 22" Back Height No Charge
 TSBH23 23" Back Height No Charge
 TSBH24 24" Back Height No Charge

CONTOURA BACK HEIGHT

- Clear
 TSBH22 22" Back Height No Charge
 TSBH25 25" Back Height No Charge
 TSBH28 28" Back Height No Charge

Back Height Chart

Back Type	Back Height
CTS - Contoura	22", 25" and 28"
BPO - Back Posts Only	20" - 24"
CURV - Curved	16" - 26"

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula® CG with MK6i™ for TDX® Spree/SC Retrofit Kit

ARM HEIGHT

ADJUSTABLE ARM TYPE (PAIR) Clear

- 19F Full Length Adj. Hgt No Charge
- 29F Desk Length Adj. Hgt No Charge

FLIP UP MIXED ADJUSTABLE ARM TYPE (EACH)

All arms are detachable.

- Right Arm** Clear
- 19FR Full Length Adj. Hgt. Right No Charge
 - 29FR Desk Length Adj. Hgt. Right No Charge

- Left Arm** Clear
- 19FL Full Length Adj. Hgt. Left No Charge
 - 29FL Desk Length Adj. Hgt. Left No Charge

CANTILEVER MIXED ARM TYPE (EACH)

- Right Arm** Clear
- CLVFR Full Length Arm Right \$175.00
 - CLVDR Desk Length Arm Right \$175.00

- Left Arm** Clear
- CLVFL Full Length Arm Left \$175.00
 - CLVDL Desk Length Arm Left \$175.00

ARMPAD OPTIONS (CANTILEVER ARMS) Clear

- WF Waterfall Armpads Standard
- RCL Reclining Armpads No Charge

ARM HEIGHT

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS Clear

- 7311BK Airline Buckle Chest Positioning Strap \$100.00 E0960
- 9515 Push Button Style Chest Positioning Strap \$100.00 E0960
- 7321BK Hook and Loop Chest Positioning Strap \$100.00 E0960

POSITIONING OPTIONS

HEADREST TYPE (OPTIONAL) Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest \$210.00 E0955
- TAR27700HW Removable & Adjustable Hardware \$210.00 E1028
- TAR128000 Two-Step Headrest \$210.00 E0955
- TAR128000HW Removable & Adjustable Hardware \$210.00 E1028
- T0022400 Large Headrest \$210.00 E0955
- T0022400HW Removable & Adjustable Hardware \$210.00 E1028
- T0022600 Small Headrest \$210.00 E0955
- T0022600HW Removable & Adjustable Hardware \$210.00 E1028

STEALTH HEADRESTS Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CP175 14" Multi Axis Headrest \$225.00 E0955
- CP175HW Removable & Adjustable Hardware \$225.00 E1028
- CP275 10" Multi Axis Headrest \$225.00 E0955
- CP275HW Removable & Adjustable Hardware \$225.00 E1028
- CP180 14" TWB Headrest \$225.00 E0955
- CP180HW Removable & Adjustable Hardware \$225.00 E1028
- CP280 10" TWB Headrest \$225.00 E0955
- CP280HW Removable & Adjustable Hardware \$225.00 E1028

LATERALS Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TCSLM Clamp-On Swingaway Laterals Medium \$210.00 E0956
- TCSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TCSLL Clamp-On Swingaway Laterals Large \$210.00 E0956
- TCSLLHW Removable and Adjustable Hardware \$210.00 E1028

STEALTH LATERALS Clear

Hardware is packaged with equipment and cannot be ordered separately.

- TWB-CPR Adj. Swing Contoura Lat. Med \$320.00 E0956
- CPRHW Removable & Adjustable Hardware \$320.00 E1028
- TWB-CPRL Adj. Swing Contoura Lat. Lg \$320.00 E0956
- CPRLHW Removable & Adjustable Hardware \$320.00 E1028

LATERALS FOR CONTOURA SEATING Clear

Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00 E0956
- CSLMHW Removable and Adjustable Hardware \$210.00 E1028
- CSSL Swingaway Contoura Lat. Lg \$210.00 E0956
- CSSLHW Removable and Adjustable Hardware \$210.00 E1028

MK6i™ ELECTRONICS

MK6i CONTROLLERS

- MK660ACC 60 Amp Controller Standard
- GTRAC G-Trac Module \$995.00 K0108

STANDARD DRIVER CONTROL

- SPJAP SPJ+ for Power Seating No Charge
- Includes egg switch with 611 mount to operate single actuator.

ACTUATOR CONTROL THROUGH SPJ+

- SPJAPM SPJ+ with Joystick Actuator Control \$1,300.00 E2310
- Allows Actuator Operation through the driver control.
Includes port for Attendant Switch.

EXPANDABLE CONTROLLER SYSTEM Clear

- PWH Harness req. for Expandable Electronics \$350.00 E2313
- Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPI+ Multiple Drive Joystick with Color Display .. \$750.00 E2377
- ASLRDYM6 MK6 Display Only (No Driver Control) \$750.00 E2377
- Display only, No Driver Control. Allows choice from order form.
- MCD Mini Color Display \$750.00 E2377
- No port for remote On/Off switch. Stereo port only.

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare® Formula® CG with MK6i™ for TDX® Spree/SC Retrofit Kit

FOR YOUR ASL DRIVE SYSTEM NEEDS
PLEASE REFER TO THE FOLLOWING:
ASL SINGLE INVOICE PRICE LIST FORM # 02-097.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS Clear

- 1500M6 Prop RIM Head Control with MK6i Display .. \$2,875.00 E2327 Clear
- SMHD Small Headrest No Charge
- LGHD Large Headrest No Charge
- 1558MM6 Compact Joystick with Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS Clear

- The below selections must also use either a Multiple Drive Joystick or MK6i Display.
- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325
 - 5020M6 Mini Tash Joystick and Digital Interface Box \$1,750.00 E2321

SIP-N-PUFF KIT

- PKG32666 Therafin® Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- 5020 Mini Tash Joystick \$405.00
Requires SNPM6 Sip-N-Puff Interface Box.
- PACM6 Prop Attendant Control with quick release mount \$850.00
- 1552M6 Digital Attendant Control \$850.00
Requires SNPM6 Sip-N-Puff Interface Box.

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ \$279.00 E1028
Include ARM250 standard.
- QLA Quad Link Retract. Joystick Mnt SPJ+ \$279.00 E1028
Include ARM250 standard.
- ARM250 Stealth Height Adj Joystick \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge
- RIMHW RIM Hardware For 1500 \$219.00 E1028
- PKG32669 Complete Bib Assembly (Tash/1558) \$265.00
- GATMPJM6 Gatlin Midline for CMPJ+ \$600.00 E1028
- GAT1812 Gatlin Midline for 1812 \$600.00 E1028
- RCM Rear Cane Mount for Compact Joystick \$395.00

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USBReady) \$90.00
- 1560 "T" Handle Flexible Joystick Ext. \$77.00 E2323
- 1561 Straight Handle Flexible Joystick Ext. \$77.00 E2323
- 1826 Chin Cup \$48.00 E2324
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for outputs 1 and 2 \$795.00
- AUX34M6 Auxiliary Module for outputs 3 and 4 \$795.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- 2SW Dual Toggle \$172.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+ Joystick options.
- PBMT Push Button Mounting Kit \$252.00
- RFM Wireless (RF) Mouse Only Module \$750.00
RFM will be available in September, 2010. Allows wireless computer mouse operation through the driver control. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ joysticks.
- IR Infrared Module with Wireless RF Mouse Control \$995.00
Capable of operating up to six remote Infra-Red controlled devices. Also allows wireless (RF) Mouse operation through the driver control. Compatible with RFM module for control of two wireless mice from one chair. Includes two button stereo switch (with hook/loop mount) for mouse clicks. Not compatible with SPJ+ Joysticks.

POWERED SEATING ACCESSORIES Clear

- SANODE One Actuator Operation thru Driver Cntrl..... \$1,300.00 E2310
Required for Tilt Operation through Driver Control - Single Actuator Systems Only. Not available with SPJ+ Joystick. Does not come with Egg switch.
- 4WSB Four Way Switch Box No Charge Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311 Clear
- FWT 4 Way Toggle No Charge
- QPB 4 Push Buttons No Charge

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Switch - Black \$156.00
- ASL304A Light Wobble Switch \$180.00
- ASL202J Fiber Optic Switch \$732.00
Not available with SPJ+ joystick. Must add MK6 Auxiliary Power Source.
- ASL208J Adjustable Proximity Switch \$660.00
Not available with SPJ+ joystick. Must add MK6 Auxiliary Power Source.
- ASL504M6 Remote Emergency Stop Switch \$1,185.00
Not available with SPJ+ joystick. Must add MK6 Auxiliary Power Source.
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS ACCESSORIES

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00

BATTERY CHARGER OPTION Clear

- 110CHARGER 110 Volt Battery Charger No Charge
- LESTR Lester 8 AMP Battery Charger \$320.00 E2367
- L8900 Omit Battery Charger Less \$80.00

VENT

- COMP Compact Vent Tray \$1,540.00
For ventilators with self-contained power supply and back up.

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

FRONT RIGGING OPTIONS

CENTER MOUNT

Clear

For a 13.25" - 17.25" extension range, available through special #SM-071217095848 with Tall STF Only and 5 Degree Seat-Tilt Angle.

CTMT-F Center Mount Adj. Knee Angle (12.25" - 15") \$300.00
 Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.

Clear

- SML Small Footplates (5.25" x 7.5") No Charge
- MDM Medium Footplates (5.5" x 9.5") No Charge
- LRG Large Footplates (5.90" x 11.5") No Charge
- CALF Calf Pads For Center Mount \$200.00

CTMT-FJ Jr. Center Mount Adj. Knee Angle (9.25" - 12") \$300.00
 Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.

Clear

- SML Small Footplates (5.25" x 7.5") No Charge
- MDM Medium Footplates (5.5" x 9.5") No Charge
- LRG Large Footplates (5.90" x 11.5") No Charge
- CALF Calf Pads For Center Mount \$200.00

CTMT-F Center Mount Extension Matrix			
0 Degree Seat Angle	Low	Med	Tall
CTMT-F	12 3/4"	13 1/2"	14 1/4"
5 Degree Seat Angle	Low	Med	Tall
CTMT-F	N/A	14 1/4"	15"

*5 degree seat angle not available with LOW seat-to-floor.

PIN STYLE - HEMI (PAIR)

Clear

Provides up to 1 1/2" of depth adjustment in 1/2" increments. PHW93 standard with 1651 composite footplate. PHAL4A standard with Aluminum footplate.

- PW93 Swingaway Footrest 13 1/2" - 19" No Charge
- PAL4A Smartleg Art Legrests 17 1/4" - 21 1/2" \$320.00 K0053

FOOTPLATE OPTIONS (PAIR)

Clear

- 1615 Composite Footplates No Charge
- 1350 Extra Lrg Aluminum Footplates No Charge
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- HLOOP Heel Loops (pair) \$45.00 E0951
- LPT2 Longer Pivot Slide Tube for 1350 \$83.00
- IMPACTG Impact Guards \$51.00
- ALPT Longer pivot And Slide Tube \$84.00
- ALPT4 Longer Pivot And Slide Tube 4" \$84.00

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula® CG with MK6i™ for TDX® Spree/SC Retrofit Kit

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Tech Tips

Operating Powered Seating through the Driver's Control. All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) - E2310.

Single Actuator Systems with CMPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310.

Features and Specifications

- Weight Capacity of 300 lbs.
- Tilt range from 0 degrees to 55 degrees or 5 degrees to 60 degrees
- Curved or Contoura style seating
- Adjustable seat frame

After Market Backrests

The following companies also have hardware to make their backrests work with our new canes:

- AEL - planar seating system
- Freedom Designs® - planar seating system
- AES - planar seating system
- Canyon - planar seating system

Contoura Seat Sizes

	D16	D17	D18	D19	D20
W16	X	X	X	X	X
W18	X	X	X	X	X
W20	X	X	X	X	X
Back Heights	22", 25", or 28"				

BPO or Curved Seat Sizes

Seat Widths	16" - 20"
Seat Depths	16" - 20"
Back Heights	20" - 26" (BPO); 16" - 26" (CURV)

Comments

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation
 USA

One Invacare Way
 Elyria, OH USA
 44035-4190
 800-333-6900

www.invacare.com

©2010 Invacare Corporation
 All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Products designated by ASL are Manufactured by Adaptive Switch Labs, Inc. Therafin is a registered trademark of W/R Medical Electronics Co. Specifications and prices are subject to change without notification.

