



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.

Please consider the environment before printing this page with your order form.
For more information on features and benefits of electronic order forms, please visit www.invacare.com/eForms

Invacare® Formula™ CG Tilt/Recline and Tilt/Recline/Elevate for TDX® SP

Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

POWERED SEATING BASE

TDXSP-MCG TDXSP for Formula CG with MK6i™ \$6,275.00 ^{K0861} _{K0862}

Group 3 - K0861 Basic Package	
TDXSP-MCG TDXSP for Formula CG with MK6i	
300 lbs. Weight Capacity	Standard
4-Pole Motors	Standard
1431-3 14" x 3" Tires with Foam Filled Inserts	No Charge
Flip Back Adjustable Height Arm	Standard
MK690ACC 90 Amp Controller	Standard
70NHD 70° Swingaway Footrests	Standard
8 Amp Off Board Battery Charger	Standard
1515 Push Button Style Seat Positioning Strap	Standard
Seat Option: CGTR Tilt/Recline (E1007)	\$9,595.00
Seat Option: CGTRE Tilt/Recline/Elevate (E1007)	\$9,595.00
Seat Option: E2-TDX Elevate to assist w/transfers (E2300)	\$1,895.00

OMIT SEAT FRAME

OSF Omit Seat Frame for Base Only Less \$700.00

TRANSPORT TIE DOWN

TRBKTS Wheelchair Transport Brackets \$250.00

TRANSPORT TIE-DOWN*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

DRIVE WHEEL OPTIONS

1430-3 14" x 3" Pneumatic Tires Standard
 1431-3 14" x 3" Tires w/Foam Filled Inserts No Charge

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

FRAME FINISH

CONTEMPORARY

- 24P Wet Black No Charge
- 60P Silver Metallic No Charge
- 105P Midnight Blue No Charge
- 115P Black Prism No Charge

TOO HOT

- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge

WAY COOL

- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 29P White Snow No Charge

SEAT

SYSTEM TYPE

- CGTR CG Tilt/Recline System with Mechanical Shear Reduction .. \$9,595.00 E1007
 - CNVACTS Conventional Actuator Standard
- CGTRE CG Tilt/Recline/Elevate \$9,595.00 E1007
 - Must choose the E2-TDX option with CGRTRE**
 - E2-TDX Elevate to assist with transfers \$1,895.00 E2300
 - CNVACTS Conventional Actuator Standard

SEAT-TO-FLOOR HEIGHT

Nominal measurement at front edge of seat pan to floor on 18" deep seat. STF bracket allows for 1" of adjustability.

- LOW 17.25" to 18.25" at 0° No Charge
- MED 18.0" to 19.0" at 0° No Charge
- TALL 19.0" to 20.0" at 0° No Charge

SEAT SELECTION

- CNB Conventional Back No Charge
- CTB Contoura™ Seating \$605.00 E2620 E2621
 - Does not include seat cushion.**

MATRIX ELITE BACK 16"-20"w (3" contour)

- PBE14 Matrix Elite 14" high \$625.00 E2620
- PBE16 Matrix Elite 16" high \$625.00 E2620
- PBE 18 Matrix Elite 18" high \$625.00 E2620
- PBE 20 Matrix Elite 20" high \$625.00 E2620

ADJUSTABLE SEAT WIDTH RANGE

- 1620W 16" - 20" Seat Width Range No Charge
- 2024W 20" - 22" Seat Width Range No Charge

SEAT WIDTH SETTING

- W16 Seat Width 16" No Charge
- W17 Seat Width 17" (Not with Contoura) No Charge
- W18 Seat Width 18" No Charge
- W19 Seat Width 19" (Not with Contoura) No Charge
- W20 Seat Width 20" No Charge
- W21 Seat Width 21" (Not with Contoura) No Charge
- W22 Seat Width 22" No Charge

ADJUSTABLE SEAT DEPTH RANGE

- 1619D 16" - 19" Seat Depth Range (Adjustable) No Charge
- 1922D 19" - 22" Seat Depth Range (Adjustable) No Charge

SEAT DEPTH SETTING

Contoura available in Seat Depths from 16" - 22".

- D16 16" Seat Depth No Charge
- D17 17" Seat Depth No Charge
- D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge
- D20 20" Seat Depth No Charge
- D21 21" Seat Depth No Charge
- D22 22" Seat Depth No Charge

SEAT TILT

- 0SA 0° Seat Tilt No Charge
- 5SA 5° Seat Tilt No Charge

MIXED ADJUSTABLE ARM TYPE (EACH)

Right Arm

- RA19FR Full Length Adj. Hgt. Right No Charge
- RA29FR Desk Length Adj. Hgt. Right No Charge

Left Arm

- RA19FL Full Length Adj. Hgt. Left No Charge
- RA29FL Desk Length Adj. Hgt. Left No Charge

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap Standard

CHEST POSITIONING STRAPS

- 7321BK Hook and Loop Chest Positioning Strap \$100.00 E0960

BACK HEIGHT

- TSBH20 20" Back Height No Charge
- TSBH21 21" Back Height No Charge
- TSBH22 22" Back Height No Charge
- TSBH23 23" Back Height No Charge
- TSBH24 24" Back Height No Charge
- TSBH25 25" Back Height No Charge
- TSBH26 26" Back Height No Charge

CONTOURA BACK HEIGHT

- TSBH25 25" Back Height No Charge
- TSBH28 28" Back Height No Charge

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

POSITIONING ACCESSORIES

HEADREST TYPE (OPTIONAL)

Hardware is packaged with equipment and cannot be ordered separately.

- TAR127700 Curved Headrest \$210.00 E0955
- TAR127700HW Removable & Adjustable Hardware \$210.00 E1028

MATRIX HEADREST PADS WITH STANDARD COVER

- MEHSP-6 Matrix Elan Standard Pad 6" \$235.00 E0955
- MEHSP-10 Matrix Elan Standard Pad 10" \$235.00 E0955
- MEHSP-14 Matrix Elan Standard Pad 14" \$235.00 E0955

MATRIX HEADREST MOUNTING HARDWARE

- MEMAH Elan Multi-Axis Removable HDW \$220.00 E1028

LATERAL HIP SUPPORTS

HIP SUPPORT & PAD - LEFT

- LHSFL Hip Support with fixed HDW Left \$120.00 E0956
- LHS-XS Size XS (5.5"Wx4"H) No Charge
- LHS-S Size S (7"Wx6"H) No Charge
- LHS-M Size M (9"Wx4"H) No Charge
- LHS-L Size L (13"Wx4"H) No Charge

HIP SUPPORT & PAD - RIGHT

- LHSFR Hip Support with fixed HDW Right \$120.00 E0956
- RHS-XS Size XS (5.5"Wx4"H) No Charge
- RHS-S Size S (7"Wx6"H) No Charge
- RHS-M Size M (9"Wx4"H) No Charge
- RHS-L Size L (13"Wx4"H) No Charge

REMOVABLE HIP SUPPORT HARDWARE

Maxx Style Quick Release, Removable, Multi-Axis Mounting HDW for Hip Supports

- MQRHL MAXX Style Quick Release - Left \$215.00 E1028
- MQRHR MAXX Style Quick Release - Right \$215.00 E1028

ADDUCTORS

Hardware is packaged with equipment and cannot be ordered separately.

- ADD Adductors Small \$200.00 E0956
- ADDHW Hardware for ADD \$400.00 E1028

LATERALS

LATERALS - CONVENTIONAL BACK

For use with Conventional seating only. Hardware is packaged with equipment and cannot be ordered separately.

- TMSLM Multi-Bar Swingaway Laterals Medium \$210.00 E0956
 - TMSLMHW Removable and Adjustable Hardware \$210.00 E1028
- TMSLL Multi-Bar Swingaway Laterals Large \$210.00 E0956
 - TMSLLHW Removable and Adjustable Hardware \$210.00 E1028

LATERALS - CONTOURA BACK

For use with Contoura seating only. Hardware is packaged with equipment and cannot be ordered separately.

- CSLM Swingaway Contoura Lat. Med \$210.00 E0956
 - CSLMHW Removable and Adjustable Hardware \$210.00 E1028
- CSLL Swingaway Contoura Lat. Lg \$210.00 E0956
 - CSLLHW Removable and Adjustable Hardware \$210.00 E1028

MATRIX FIXED LATERALS (FOR ELITE BACKS)

- PBLS-L Matrix Std Lateral - Fixed - Left (7.5"x5") \$100.00 E0956
- PBLS-R Matrix Std. Lateral - Fixed - Right (7.5"x5") .. \$100.00 E0956
- PBOL-L Matrix Offset Lateral Fixed Left \$110.00 E0956
- PBOL-R Matrix Offset Lateral Fixed Right \$110.00 E0956

MATRIX SWING AWAY LATERAL PADS (FOR ELITE BACKS)

LEFT TRUNK SUPPORT PAD

- L-XS Size XS (3.25"Wx4.25"H) \$100.00 E0956
- L-S Size S (3.25"Wx5.25"H) \$100.00 E0956
- L-M Size M (3.75"Wx5.75"H) \$100.00 E0956
- L-L Size L (4.25"Wx6.5"H) \$100.00 E0956

RIGHT TRUNK SUPPORT PAD

- R-XS Size XS (3.25"Wx4.25"H) \$100.00 E0956
- R-S Size S (3.25"Wx5.25"H) \$100.00 E0956
- R-M Size M (3.75"Wx5.75"H) \$100.00 E0956
- R-L Size L (4.25"Wx6.5"H) \$100.00 E0956

SWING AWAY HARDWARE (FOR ELITE BACKS)

- SAEHBL Matrix Elite Swing Away HDW Left \$195.00 E1028
- SAHEBR Matrix Elite Swing Away HDW Right \$195.00 E1028
- OSAHL Matrix Offset Elite Swing Away HDW Left ... \$195.00 E1028
- OSAHR Matrix Offset Elite Swing Away HDW Right . \$195.00 E1028

NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.

ELECTRONICS

MK6i CONTROLLER

- MK690ACC Controller with Actuator Control Standard
- GTRAC G-Trac Module \$995.00 K0108
G-Trac provides improved tracking capabilities and enhanced traction. Does not require special motors.

EXPANDABLE CONTROLLER SYSTEM

- PWH Harness req. for Expandable Electronics \$350.00 E2313
Required to interface joystick or MK6i display with expandable controllers.
- CMPJM6U MPJ+™ Multiple Drive Joystick with Color Display \$750.00 E2377
- ASLRDYM6 MK6i Display Only (No Driver Control) \$750.00 E2377
Display only, No Driver Control. Allows choice from order form.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- 1558MM6 Compact Joystick w/Display \$1,130.00 E2373

ALTERNATIVE DIGITAL CONTROL SYSTEMS

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Digital Interface Box \$1,500.00 E2325

SIP-N-PUFF KIT

- PKG32666 Therafin Sip-N-Puff Breath Tube Kit \$380.00 E2326

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Std Mount \$770.00
- PACM6 Prop Attendant Control with quick release mount \$850.00

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retract. Joystick Mnt CMPJ+ \$279.00 E1028
QLAM6 includes ARM250 standard.
- ARM250 Stealth Height Adjustable Joystick Mount \$210.00 E1028
- LEFT Left Handed Mounted Joystick No Charge

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Pro Memory Card \$475.00
- 1813SD Pro Memory Card (USBReady) \$90.00
- A24VPS 24V MK6 Auxiliary Power Source \$200.00
- AUX12M6 Auxiliary Module for 1 and 2 \$795.00
- AUX34M6 Auxiliary Module for 3 and 4 \$795.00
- MPY Monoport "Y" Cable or CMPJ+ or MK6 Display .. \$15.00
- 2PB Dual Push Button \$147.00
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Not compatible with SPJ+™ Joystick options.
- PBMT Push Button Mounting Kit \$252.00

POWERED SEATING ACCESSORIES

- 4WSB Four Way Switch Box No Charge
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge
- S4WSB Multiple Actuator Interface Box \$2,600.00 E2311
 - FWT 4 Way Toggle No Charge
 - QPB 4 Push Buttons No Charge

Systems

For all systems with two or more actuators you must choose 4WSB at No Charge or the S4WSB with and Up - Charge.

S4WSB Replaces the Four Way Switch Box and allows multiple actuator operations through the driver control.

FWT or QPB is required when choosing the 4WSB OR S4WSB.

ASL OPTIONAL ACCESSORIES

- EGSBLK Egg Reset-Black \$156.00
- ASL611 Adjustable Switch Mounting Bracket \$252.00

MISCELLANEOUS

BATTERY TRAY TYPE

22TRY Base Width is 24"; 24TRY Base Width is 25 1/2".

- 22TRY 22NF Style Trays Standard
- 24TRY 24 Group Style Tray No Charge

BATTERIES

- 22NFBATTERY 22NF Battery (On Chair) \$300.00 E2361
- 24BATTERY 24Gel Battery (On Chair) \$410.00 E2363

BATTERY CHARGER OPTION

- 110CHARGER 110 Volt Battery Charger No Charge
- L8900 Omit Battery Charger Less \$80.00

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$45.00

BACK ACCESSORIES

- L8601 Omit Back Upholstery No Charge
- RBP Reinforced Back Upholstery No Charge

**NOTE: All specifications and dimensions are approximate.
 For the most current pricing information, go to www.invacare.com/homecare.**

FRONT RIGGINGS

HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket No Charge

CENTER PIVOT STYLE (PAIR)

- 70NHD 70° Swingaway Footrests 13 1/2" - 17 1/2" Standard
- 70STAPER 70° Tapered Footrests 13" - 17" No Charge
Standard with impact guards. Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.
- AT5044 Manual Elevating Legrests \$255.00 E0990

CENTER MOUNT

- CTMT-F Center Mount Adj. Knee Angle (12 1/4" - 15") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 12 1/2" - 15". Footplates are adjustable angle. Must choose either SM, MED or LRG footplate.
 - SM Small Footplates (5 1/4" x 7 1/2") No Charge
 - MED Medium Footplates (5 1/2" x 9 1/2") No Charge
 - LRG Large Footplates (5.90" x 11 1/2") No Charge
 - CALF Calf Pads For Center Mount \$200.00
- CTMT-FJ Jr. Center Mount Adj. Knee Angle (9 1/4" - 12") \$300.00
Allows 60° to 90° static knee angle positioning with a drop back lever for easier transferring. Extension range 9 1/4" - 12". Footplates are adjustable angle. Available on low seat-to-floor only. Must choose either SM, MED or LRG footplate.
 - SM Small Footplates (5 1/4" x 7 1/2") No Charge
 - MED Medium Footplates (5 1/2" x 9 1/2") No Charge
 - LRG Large Footplates (5.90" x 11 1/2") No Charge
 - CALF Calf Pads For Center Mount \$200.00

POWER FOOT POSITIONING OPTIONS

- LNX and LNX-L come Standard with curved individual calf pads that are independently height and angle adjustable. KtH range: LNX=13" to 15.5", LNX-L=13" to 19"
- LNX Power Center Mount Foot Platform \$2,995.00 E1012
Comes standard with a footboard. Not available with SPJ+ Joystick.
 - SML Small Footplates (5.25" x 7.5") \$275.00
SML comes with 2" of width adjustment ranging from 12-14"
 - MDM Medium Footplates (5.5" x 9.5") \$275.00
MDM comes with 2" of width adjustment ranging from 13-15"
 - LRG Large Footplates (5.9" x 11.5") \$275.00
LRG comes with 2" of width adjustment ranging from 14-16"
 - LNX-L Power Center Mount Foot Platform- Long \$2,995.00 E1012
If OSA is selected, LNX-L is only available with Tall STF height. Comes standard with a footboard. Not available with SPJ+ Joystick.
 - SML Small Footplates (5.25" x 7.5") \$275.00
SML comes with 2" of width adjustment ranging from 12-14"
 - MDM Medium Footplates (5.5" x 9.5") \$275.00
MDM comes with 2" of width adjustment ranging from 13-15"
 - LRG Large Footplates (5.9" x 11.5") \$275.00
LRG comes with 2" of width adjustment ranging from 14-16"

FOOTPLATE OPTIONS (PAIR)

- 1651 Flip-Up Composite Footplates No Charge
- 1350 Extra Large Aluminum Footplates No Charge
- AT5543 Adjustable Angle Flip-Up Footplates \$165.00 K0040

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap \$27.00 E0038
- IMPACTG Impact Guards \$51.00
- HLOOP Heel Loops (pair) \$45.00 E0951
- LPT2 Longer Pivot Slide Tube For 1350 \$83.00
- ALPT Longer Pivot and Slide Tube \$84.00
- ALPT4 Longer Pivot and Slide Tube 4" \$84.00

Contoura Seat Sizes							
	D16	D17	D18	D19	D20	D21	D22
W16	*	*	*	*	*	*	*
W18	*	*	*	*	*	*	*
W20	*	*	*	*	*	*	*
W22	*	*	*	*	*	*	*
*Not available with Contoura Seat Cushion							
Back Heights	22", 25" or 28"						

Conventional Seat Sizes	
Seat Widths	16" - 22"
Seat Depths	16" - 22"
Back Heights	20" - 26"

**NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.**

Invacare® Formula™ CG Tilt/Recline and Tilt/Recline/Elevate for TDX® SP Power Wheelchair

Tech Tips

Operating Powered Seating through the Driver's Control.
All chairs require adding or replacing one component to allow actuator operation through the driver control.

SPJ+ with PSS:

Replace joystick with SPJ+ with ACC (Actuator Control) - E2310.

Single Actuator Systems with CMPJ+, PSR+, PSF+ MK6i Display:

Add a SANODE (Single Actuator Node) to the system - E2310.

Multiple Actuator Systems:

Replace the standard 4 Way Switch box with a S4WSB (Multiple Actuator Interface Box) - E2311.

Formula CG Tilt/Recline Features and Specifications

- Weight capacity of 300 lbs. or 400 lbs. (with heavy duty 4-pole motors).
- Conventional or Contoura style seating.
- Adjustable seat frame with Conventional style.

After Market Backrests

The following companies also have hardware to make their backrests work with our new canes:

- AEL - planar seating system
- Freedom Designs® - planar seating system
- Canyon - planar seating system

Seating and Positioning

For your seating needs, please refer to the following rehab seating and positioning products: Price List Form No. 07-027

Batteries (not included)

22NF batteries with terminal configuration (positive on the left and negative on the right) **MUST** be used. 22NF batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) **MUST** be used. GP24 batteries that have the reverse terminal configuration **MUST** not be used. Terminals **MUST** have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations **MUST** be followed otherwise injury and damage may occur.

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, PSR+, PSF+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

Comments

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way

Elyria, OH USA

44035-4190

800-333-6900

©2018 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.



Yes, you can.®

Revised July 2, 2018