#### **INVACARE®**

# **PINDOT**

**CUSTOM SEATING** 



## **Personalized Seating Solutions**

State-of-the-art technology and computer-assisted manufacturing from PinDot produces custom contoured seat and back cushions for precise support and excellent stability when off-the-shelf products wont accomplish your seating needs.

Call to order **1.800.451.3553**invacare.com







## PINDOT® CONTOURU® & CONTOURU® 4U

With PinDot ContourU, user simulation creates an individually shaped support surface for optimal support and stability and a precise individual fit that provides exceptional pressure redistribution and positioning.

Multiple finishing options are available, including "naked" foam with a removable cover for increased comfort or a non-removable vacuformed vinyl for stability and increased foam protection.

ContourU 4U is a pre-configured solution for ease of ordering. Simply send your 3D CLIP image, pick your mounting, width, everything else is pre-selected for you. Quicker quotes, orders and delivery.

Aluminum or ABS plastic reinforcements are custom fit and built into the cushion, minimizing the need for additional support accessories.



## PINDOT® SILHOUETTE® SERIES

PinDot Silhouette creates optimal product contouring for positioning and balance within the wheelchair frame through User Simulation or Molding. Precise, custom-carved polyurethane foam in soft, firm or extra firm provides exceptional pressure redistribution and positioning, even for users with orthopedic asymmetries.


A moisture-proof coating seals and protects the foam from moisture and allows the cushions to be easily wiped clean. For additional foam protection, add an optional moisture-proof and anti-microbial inner cover.


Utilize the Invacare® Matrx® Elite™, Elite™ TR and Elite™ PB back shells for a more contoured look with the ease and adjustability of the Matrx Easy Set Hardware.




### PINDOT® SILHOUETTE CUSTOM BASIC


The PinDot Silhouette Custom Basic sets simulation aside and allows you to design a cushion to fit the needs of your client. Start with pre-captured shapes, use measurements from your mat evaluation and specify any needed modifications through our computer-aided design software, or fill out an order form and let PinDot do the work for you.












**Pre-captured Shape** 





To request a demo visit invacare.com/clip or call 1.800.451.3553 elyriacustomseating@invacare.com



















**Invacare Corporation** www.invacare.com

1200 Taylor Street Elyria, Ohio 44035-4190 (800) 333-6900

#### Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

#### ©2018

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 18-034