

FREEDOM® DESIGNS
STROLLER HANDLE
REMOVABLE & ADJUSTABLE



Freedom® Designs is reintroducing this very popular and sought after wheelchair accessory, the **Removable & Adjustable Stroller Handle** with improved features and enhanced functionality.

Call to order

1.800-331-8551
Freedomdesigns.com



Yes, you can.®



Our kids are our passion.™

Freedom Designs Removable & Adjustable Stroller Handle Features and Functional Benefits

This bend provides enough height to push the chair while in a moderate amount of posterior tilt. While the hand positioning creates a comfortable distance from the rear of the wheelchair making room for ADLs and other accessory items.

Developed with ease of use in mind, the handle features a slight amount of lateral movement built-in to make installation and removal on a tilt-in-space wheelchair frame quick and easy.

The 115 degree bend places the caregivers hands at a comfortable ergonomic angle.



- Width-Adjustable (2" Adjustment, Grows with frame)

- Includes black neoprene wrap pad

- Removable, detent-pin release

Tilt-in-space wheelchair sold separately

- Reimbursement friendly alternative
- Growable center section creates a flex-point to make installation and removal a breeze on tilt-in-space frames
- Back post tube inside diameter is reamed .031" for ease of installation, requiring no additional bushings for fitment
- Retrofittable to existing chairs with service parts
- Designed for the Freedom® NXT® and is compatible with all Freedom Designs frames, including the Freedom® P.R.O. CG™ Tilt-in-Space Wheelchair
 - *Not compatible with recline back on Freedom® P.R.O. CG™ Tilt-in-Space Wheelchair
- Compatible with standard 90° bend or Ergo bend push handles

Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Freedom Designs
www.freedomdesigns.com

USA
2241 N. Madera Rd.
Simi Valley, CA 93065
(800) 331-8551
Fax: (888) 582-1509

© 2018
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
©. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 18-273
180929

